

JOBST[®] Ready-to-Wear

Ordering Information



JOBST[®] Bella[™] Lite
JOBST[®] Bellavar[®]
JOBST[®] Opaque
JOBST[®] UltraSheer

Contents

Ready-to-Wear Compression Garments

Order Codes





JOBST® Bella™ Lite	3
JOBST® Bellavar®	4-5
JOBST® Opaque	6-25
JOBST® UltraSheer	26-39

Measuring Guides

JOBST® Bella™ Lite	3
JOBST® Bellavar®	5
JOBST® Opaque	39
JOBST® UltraSheer	39
Contact us	40

Please note: we are unable to include all styles, options and colours for certain products therefore if you cannot find the exact product you require, please contact our Customer Service team for further assistance. Additional information can also be found on our website at www.bsnmedical.co.uk

Ordering Information

	BSN Code	Option	Size	Length	Colour	PIP	NHSSC	
Armsleeve CCL1 	75273-02	Knitted Band	Small	Regular	Beige	382-1196	EGD7879	
	75273-00	Knitted Band	Small	Long	Beige	325-7011	EGD340	
	75274-02	Knitted Band	Medium	Regular	Beige	382-1188	EGD7880	
	75274-00	Knitted Band	Medium	Long	Beige	325-7037	EGD341	
	75275-02	Knitted Band	Large	Regular	Beige	382-1170	EGD7881	
	75275-00	Knitted Band	Large	Long	Beige	325-7029	EGD342	
	75273-03	Silicone Band	Small	Regular	Beige	382-1238	EGD7882	
	75273-01	Silicone Band	Small	Long	Beige	325-6955	EGD346	
	75274-03	Silicone Band	Medium	Regular	Beige	382-1220	EGD7883	
	75274-01	Silicone Band	Medium	Long	Beige	325-6971	EGD347	
	75275-03	Silicone Band	Large	Regular	Beige	382-1212	EGD7884	
	75275-01	Silicone Band	Large	Long	Beige	325-6963	EGD348	
Armsleeve CCL2 	75276-05	Knitted Band	Small	Regular	Beige	382-1261	EGD7885	
	75276-00	Knitted Band	Small	Long	Beige	325-6989	EGD343	
	75277-05	Knitted Band	Medium	Regular	Beige	382-1204	EGD7886	
	75277-00	Knitted Band	Medium	Long	Beige	325-7003	EGD344	
	75278-05	Knitted Band	Large	Regular	Beige	382-1246	EGD7887	
	75278-00	Knitted Band	Large	Long	Beige	325-6997	EGD345	
	75276-06	Silicone Band	Small	Regular	Beige	382-1287	EGD7888	
	75276-01	Silicone Band	Small	Long	Beige	325-6898	EGD7891	
	75277-06	Silicone Band	Medium	Regular	Beige	382-1279	EGD7889	
	75277-01	Silicone Band	Medium	Long	Beige	325-6914	EGD350	
	75278-06	Silicone Band	Large	Regular	Beige	382-1253	EGD7890	
	75278-01	Silicone Band	Large	Long	Beige	325-6906	EGD351	
	75276-04	Knitted Band	Small	Long	Black	363-7568		
	75277-04	Knitted Band	Medium	Long	Black	363-7576		
	75278-04	Knitted Band	Large	Long	Black	363-7550		
	Gauntlet CCL1 	75279-00	With Thumb	Small	-	Beige	325-6823	EGD352
75280-00		With Thumb	Medium	-	Beige	325-6831	EGD353	
75281-00		With Thumb	Large	-	Beige	325-6922	EGD354	
Gauntlet CCL2 		75282-00	With Thumb	Small	-	Beige	325-6773	EGD355
		75283-00	With Thumb	Medium	-	Beige	325-6781	EGD356
		75284-00	With Thumb	Large	-	Beige	325-6799	EGD357
		75282-03	With Thumb	Small	-	Black	363-7543	
		75283-03	With Thumb	Medium	-	Black	363-7535	
		75284-03	With Thumb	Large	-	Black	363-7527	

Gauntlet

Measuring Guide

Armsleeve				Gauntlet			Length
							
Size	Wrist	Mid-Lower	Mid-Upper	Size	Palm	Wrist	Note: length measurement taken on inside of arm
Small	14cm-20cm	18cm-23cm	22cm-40cm	Small	15cm-19cm	14cm-20cm	Regular up to (≤) 43cm
Medium	16cm-22cm	23cm-29cm	27cm-45cm	Medium	19cm-22cm	16cm - 22cm	Long over (>) 43cm
Large	18cm-24cm	29cm-34cm	33cm-51cm	Large	22cm - 26cm	18cm - 24cm	

Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Knee High - Sahara ■

RAL CCL 2 - Open Toe

Petite			Regular			
I	77570-00	326-1195	EGD5028	77571-00	326-1138	EGD5034
II	77570-01	326-1203	EGD5029	77571-01	326-1146	EGD5035
III	77570-02	326-1211	EGD5030	77571-02	326-1153	EGD5036
IV	77570-03	326-1229	EGD5037	77571-03	326-1161	EGD5043
V	77570-04	326-1237	EGD5032	77571-04	326-1179	EGD5038
VI	77570-05	326-1245	EGD5033	77571-05	326-1187	EGD5039

RAL CCL 3 - Open Toe

Petite			Regular		
I	77581-12	326-1377	77582-00	326-1435	EGD5071
II	77581-13	326-1385	77582-01	326-1443	EGD4957
III	77581-14	326-1393	77582-02	326-1450	EGD4964
IV	77581-15	326-1401	77582-03	326-1468	EGD4959
V	77581-16	326-1419	77582-04	326-1476	EGD4960
VI	77581-17	326-1427	77582-05	326-1484	EGD4961

RAL CCL 2 - Closed Toe

Petite, short foot		Regular, short foot		Petite, long foot		Regular, long foot	
I	77570-12	77571-12	77570-24	77571-24			
II	77570-13	77571-13	77570-25	77571-25			
III	77570-14	77571-14	77570-26	77571-26			
IV	77570-15	77571-15	77570-27	77571-27			
V	77570-16	77571-16	77570-28	77571-28			
VI	77570-17	77571-17	77570-29	77571-29			

RAL CCL 3 - Closed Toe

Petite, short foot		Regular, short foot		Petite, long foot		Regular, long foot	
I	77581-24	77582-12	77581-36	77582-24			
II	77581-25	77582-13	77581-37	77582-25			
III	77581-26	77582-14	77581-38	77582-26			
IV	77581-27	77582-15	77581-39	77582-27			
V	77581-28	77582-16	77581-40	77582-28			
VI	77581-29	77582-17	77581-41	77582-29			

Thigh High - Sahara ■

RAL CCL 2 - Open Toe

Knitted band			Dotted silicone band		
I	77573-00	326-1252	77573-33	326-1310	EGD5012
II	77573-01	326-1260	77573-11	326-1328	EGD5013
III	77573-35	326-1278	77573-12	326-1336	EGD5008
IV	77573-02	326-1286	77573-13	326-1344	EGD5009
V	77573-03	326-1294	77573-14	326-1351	EGD5010
VI	77573-04	326-1302	77573-15	326-1369	EGD5011

RAL CCL 3 - Open Toe

Knitted band			Dotted silicone band		
I	77584-12	364-0380	77584-24	364-0448	
II	77584-13	364-0398	77584-25	364-0455	
III	77584-14	364-0406	77584-26	364-0463	
IV	77584-15	364-0414	77584-27	364-0471	
V	77584-16	364-0422	77584-28	364-0489	
VI	77584-17	364-0430	77584-29	364-0497	

Tights - Sahara ■

RAL CCL 2 - Open Toe

Petite			Regular		
I	77577-20	364-0786	77576-08	364-0729	
II	77577-21	364-0794	77576-09	364-0737	
III	77577-22	364-0802	77576-10	364-0745	
IV	77577-23	364-0810	77576-11	364-0752	
V	77577-08	364-0828	77576-12	364-0760	
VI	77577-0	364-0836	77576-13	364-0778	

RAL CCL 3 - Open Toe

Petite			Regular		
I	77588-00	364-0562	77587-00	364-0505	
II	77588-01	364-0570	77587-01	364-0513	
III	77588-02	364-0588	77587-02	364-0521	
IV	77588-03	364-0596	77587-03	364-0539	
V	77588-04	364-0604	77587-04	364-0547	
VI	77588-05	364-0612	77587-05	364-0554	

Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Knee High - Pearl ■

RAL CCL 2 - Open Toe

Petite			Regular		
I	77570-58	383-5246	77571-36	383-5188	
II	77570-59	383-5261	77571-37	383-5196	
III	77570-60	383-5279	77571-38	383-5204	
IV	77570-61	383-5287	77571-39	383-5212	
V	77570-62	383-5295	77571-40	383-5220	
VI	77570-63	383-5303	77571-41	383-5238	

RAL CCL 3 - Open Toe

Petite			Regular		
I	77581-58	383-5410	77582-36	383-5352	
II	77581-59	383-5428	77582-37	383-5360	
III	77581-60	383-5436	77582-38	383-5378	
IV	77581-61	383-5444	77582-39	383-5386	
V	77581-62	383-5451	77582-40	383-5394	
VI	77581-63	383-5469	77582-41	383-5402	

Thigh High - Pearl ■

RAL CCL 2 - Open Toe

Knitted band			Dotted silicone band		
I	77573-42	383-5550	77573-48	383-5626	
II	77573-43	383-5568	77573-49	383-5634	
III	77573-44	383-5576	77573-50	383-5642	
IV	77573-45	383-5584	77573-51	383-5659	
V	77573-46	383-5584	77573-52	383-5659	
VI	77573-47	383-5600	77573-53	383-5675	

RAL CCL 3 - Open Toe

Knitted band			Dotted silicone band		
I	77584-42	383-6319	77584-48	383-6251	
II	77584-43	383-6327	77584-49	383-6269	
III	77584-44	383-6335	77584-50	383-6277	
IV	77584-45	383-6343	77584-51	383-6285	
V	77584-46	383-6350	77584-52	383-6293	
VI	77584-47	383-6368	77584-53	383-6301	

Tights - Pearl ■

RAL CCL 2 - Open Toe

Petite			Regular		
I	77577-30	383-5998	77576-30	383-5931	
II	77577-31	383-6004	77576-31	383-5949	
III	77577-32	383-6012	77576-32	383-5956	
IV	77577-33	383-6020	77576-33	383-5964	
V	77577-34	383-6038	77576-34	383-5972	
VI	77577-35	383-6046	77576-35	383-5980	

RAL CCL 3 - Open Toe

Petite			Regular		
I	77588-24	383-6111	77587-24	383-6053	
II	77588-25	383-6129	77587-25	383-6061	
III	77588-26	383-6137	77587-26	383-6079	
IV	77588-27	383-6145	77587-27	383-6087	
V	77588-28	383-6152	77587-28	383-6095	
VI	77588-29	383-6152	77587-29	383-6103	

Measuring Guide

Measurement Point	circumference in cm						Measurement Point	circumference in cm					
	I	II	III	IV	V	VI		I	II	III	IV	V	VI
cG	42-44	46-48	50-52	54-56	58-60	62-64	cB1	23-24	25-26	27-28	29-30	31-32	33-34
cF	39-41	42-44	45-47	48-50	51-53	54-56	cB	18-19	20-21	22-23	24-25	26-27	28-29
cE	29-31	32-34	35-37	38-40	41-43	44-46	cY	27-28	29-30	31-32	33-34	35-36	37-38
cD	26-28	29-31	32-34	35-37	38-40	41-43	cA	18-19	20-21	22-23	24-25	26-27	28-29
cC	29-30	31-33	34-36	37-39	40-42	43-45							

Length measurements for all sizes (in cm):
 AD: ID petite 36-39, ID regular 40-43 AG: IG standard 69-74 AT/AG-T: IK petite 70-75, IK regular 76-80 AD: Shoe sizes with closed toe: short - up to shoe size 39 long - from shoe size 40



Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Thigh High - Sand ■

RAL CCL 1 - Open Toe, Lace Silicone Band

	Petite	Petite, wider circumference	Regular	Regular, wider circumference
I	70618-192 357-3896	70618-210 357-4316	70617-195 357-4258	70617-213 357-4373
II	70618-193 357-3904	70618-211 357-4324	70617-196 357-4266	70617-214 357-4381
III	70618-194 357-3912	70618-212 357-4332	70617-197 357-4274	70617-215 357-4399
IV	70618-195 357-3920	70618-213 357-4340	70617-198 357-4282	70617-216 357-4407
V	70618-196 357-3938	70618-214 357-4357	70617-199 357-4290	70617-217 357-4415
VI	70618-197 357-3946	70618-215 357-4365	70617-200 357-4308	70617-218 357-4423

Thigh High - Sand ■

RAL CCL 1 - Closed Toe, Lace Silicone Band

	Petite	Petite, wider circumference	Regular	Regular, wider circumference
I	70618-42 357-2518	70618-126 357-3078	70617-42 357-2658	70617-126 357-3201
II	70618-43 357-2526	70618-127 357-3094	70617-43 357-2666	70617-127 357-3219
III	70618-44 357-2542	70618-128 357-3102	70617-44 357-2674	70617-128 357-3235
IV	70618-45 357-2559	70618-129 357-3128	70617-45 357-2690	70617-129 357-3250
V	70618-46 357-2567	70618-130 357-3144	70617-46 357-2708	70617-130 357-3268
VI	70618-47 357-2575	70618-131 357-3169	70617-47 357-2724	70617-131 357-3284

Thigh High - Black ■

RAL CCL 1 - Open Toe, Dotted Silicone Band

	Petite	Petite, wider circumference	Regular	Regular, wider circumference
I	70618-246 363-8772	70618-378 363-8897	70617-249 363-8830	70617-381 363-8954
II	70618-247 363-8780	70618-379 363-8905	70617-250 363-8848	70617-382 363-8962
III	70618-248 363-8798	70618-380 363-8913	70617-251 363-8855	70617-383 363-8970
IV	70618-249 363-8806	70618-381 363-8921	70617-252 363-8863	70617-384 363-8988
V	70618-250 363-8814	70618-382 363-8939	70617-253 363-8871	70617-385 363-8996
VI	70618-251 363-8822	70618-383 363-8947	70617-254 363-8889	70617-386 363-9002

Thigh High - Black ■

RAL CCL 1 - Closed Toe, Dotted Silicone Band

	Petite	Petite, wider circumference	Regular	Regular, wider circumference
I	70618-282 363-9010	70618-414 363-9762	70617-285 363-9077	70617-417 363-9846
II	70618-283 363-9028	70618-415 363-9770	70617-286 363-9085	70617-418 363-9853
III	70618-284 363-9036	70618-416 363-9788	70617-287 363-9093	70617-419 363-9861
IV	70618-285 363-9044	70618-417 363-9796	70617-288 363-9101	70617-420 363-9879
V	70618-286 363-9051	70618-418 363-9804	70617-289 363-9119	70617-421 363-9887
VI	70618-287 363-9069	70618-419 363-9812	70617-290 363-9127	70617-422 363-9895

Thigh High - Black ■

RAL CCL 1 - Open Toe, Soft Silicone Band

	Petite	Petite, wider circumference	Regular	Regular, wider circumference
I	70618-168 327-5492	70618-174 327-5351	70617-168 327-5773	70617-174 327-5856
II	70618-169 327-5658	70618-175 327-5377	70617-169 327-5781	70617-175 327-5880
III	70618-170 327-5674	70618-176 327-5385	70617-170 327-5807	70617-176 327-5898
IV	70618-171 327-5708	70618-177 327-5393	70617-171 327-5823	70617-177 327-5914
V	70618-172 327-5724	70618-178 327-5401	70617-172 327-5831	70617-178 327-5930
VI	70618-173 327-5732	70618-179 327-5427	70617-173 327-5849	70617-179 327-5948

Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Thigh High - Black ■

RAL CCL 1 - Closed Toe, Soft Silicone Band

	Petite	Petite, wider circumference	Regular	Regular, wider circumference
I	70618-48 357-1866	70618-132 357-2146	70617-48 357-1981	70617-132 357-2237
II	70618-49 357-1874	70618-133 357-2161	70617-49 357-1999	70617-133 357-2245
III	70618-50 357-1882	70618-134 357-2179	70617-50 357-2005	70617-134 357-2260
IV	70618-51 357-1890	70618-135 357-2195	70617-51 357-2013	70617-135 357-2278
V	70618-52 357-1908	70618-136 357-2203	70617-52 357-2021	70617-136 357-2294
VI	70618-53 357-1916	70618-137 357-2211	70617-53 357-2039	70617-137 357-2310

Thigh High - Black ■

RAL CCL 1 - Open Toe, Lace Silicone Band

	Petite	Petite, wider circumference	Regular	Regular, wider circumference
I	70618-186 357-3557	70618-204 357-4084	70617-189 357-3953	70617-207 357-4027
II	70618-187 357-3581	70618-205 357-4092	70617-190 357-3961	70617-208 357-4035
III	70618-188 357-3599	70618-206 357-4100	70617-191 357-3979	70617-209 357-4043
IV	70618-189 357-3615	70618-207 357-4118	70617-192 357-3987	70617-210 357-4050
V	70618-190 357-3631	70618-208 357-4126	70617-193 357-3995	70617-211 357-4068
VI	70618-191 357-3649	70618-209 357-4134	70617-194 357-4019	70617-212 357-4076

Thigh High - Black ■

RAL CCL 1 - Closed Toe, Lace Silicone Band

	Petite	Petite, wider circumference	Regular	Regular, wider circumference
I	70618-54 357-2591	70618-138 357-3300	70617-54 357-2732	70617-138 357-3441
II	70618-55 357-2609	70618-139 357-3326	70617-55 357-2757	70617-139 357-3458
III	70618-56 357-2617	70618-140 357-3334	70617-56 357-2773	70617-140 357-3466
IV	70618-57 357-2625	70618-141 357-3342	70617-57 357-2799	70617-141 357-3474
V	70618-58 357-2633	70618-142 357-3409	70617-58 357-2815	70617-142 357-3482
VI	70618-59 357-2641	70618-143 357-3433	70617-59 357-2823	70617-143 357-3490

Thigh High - Navy ■

RAL CCL 1 - Open Toe, Dotted Silicone Band

	Petite	Petite, wider circumference	Regular	Regular, wider circumference
I	70618-240 363-9903	70618-372 364-0026	70617-243 363-9960	70617-375 364-0083
II	70618-241 363-9911	70618-373 364-0034	70617-244 363-9978	70617-376 364-0091
III	70618-242 363-9929	70618-374 364-0042	70617-245 363-9986	70617-377 364-0109
IV	70618-243 363-9937	70618-375 364-0059	70617-246 363-9994	70617-378 364-0117
V	70618-244 363-9945	70618-376 364-0067	70617-247 363-0000	70617-379 364-0125
VI	70618-245 363-9952	70618-377 364-0075	70617-248 363-0018	70617-380 364-0133

Thigh High - Navy ■

RAL CCL 1 - Closed Toe, Dotted Silicone Band

	Petite	Petite, wider circumference	Regular	Regular, wider circumference
I	70618-270 364-0141	70618-402 364-0265	70617-273 364-0208	70617-405 364-0323
II	70618-271 364-0158	70618-403 364-0273	70617-274 364-0216	70617-406 364-0331
III	70618-272 364-0166	70618-404 364-0281	70617-275 364-0224	70617-407 364-0349
IV	70618-273 364-0174	70618-405 364-0299	70617-276 364-0232	70617-408 364-0356
V	70618-274 364-0182	70618-406 364-0307	70617-277 364-0240	70617-409 364-0364
VI	70618-275 364-0190	70618-407 364-0315	70617-278 364-0257	70617-410 364-0372



Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Thigh High - Navy ■ RAL CCL 1 - Open Toe, Soft Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-216 360-8593	70618-348 360-8650	70617-219 360-8718	70617-351 360-8775
II 70618-217 360-8601	70618-349 360-8668	70617-220 360-8726	70617-352 360-8783
III 70618-218 360-8619	70618-350 360-8676	70617-221 360-8734	70617-353 360-8791
IV 70618-219 360-8627	70618-351 360-8684	70617-222 360-8742	70617-354 360-8809
V 70618-220 360-8635	70618-352 360-8692	70617-223 360-8759	70617-355 360-8817
VI 70618-221 360-8643	70618-353 360-8700	70617-224 360-8767	70617-356 360-8825

Thigh High - Navy ■ RAL CCL 1 - Closed Toe, Soft Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-24 360-7611	70618-108 360-7678	70617-24 360-7736	70617-108 360-7793
II 70618-25 360-7629	70618-109 360-7686	70617-25 360-7744	70617-109 360-7801
III 70618-26 360-7637	70618-110 360-7694	70617-26 360-7751	70617-110 360-7819
IV 70618-27 360-7645	70618-111 360-7702	70617-27 360-7769	70617-111 360-7827
V 70618-28 360-7652	70618-112 360-7710	70617-28 360-7777	70617-112 360-7835
VI 70618-29 360-7660	70618-113 360-7728	70617-29 360-7785	70617-113 360-7843

Thigh High - Navy ■ RAL CCL 1 - Open Toe, Lace Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-222 360-8833	70618-354 360-8890	70617-225 360-8957	70617-357 360-9013
II 70618-223 360-8841	70618-355 360-8908	70617-226 360-8965	70617-358 360-9021
III 70618-224 360-8858	70618-356 360-8916	70617-227 360-8973	70617-359 360-9039
IV 70618-225 360-8866	70618-357 360-8924	70617-228 360-8981	70617-360 360-9047
V 70618-226 360-8874	70618-358 360-8932	70617-229 360-8999	70617-361 360-9054
VI 70618-227 360-8882	70618-359 360-8940	70617-230 360-9005	70617-362 360-9062

Thigh High - Navy ■ RAL CCL 1 - Closed Toe, Lace Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-30 360-7850	70618-114 360-7918	70617-30 360-7975	70617-114 360-8031
II 70618-31 360-7868	70618-115 360-7926	70617-31 360-7983	70617-115 360-8049
III 70618-32 360-7876	70618-116 360-7934	70617-32 360-7991	70617-116 360-8056
IV 70618-33 360-7884	70618-117 360-7942	70617-33 360-8007	70617-117 360-8064
V 70618-34 360-7892	70618-118 360-7959	70617-34 360-8015	70617-118 360-8072
VI 70618-35 360-7900	70618-119 360-7967	70617-35 360-8023	70617-119 360-8080

Thigh High - Bronze ■ RAL CCL 1 - Open Toe, Dotted Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-444 383-8711	70618-492 383-9073	70617-447 383-7069	70617-495 383-7630
II 70618-445 383-8729	70618-493 383-9081	70617-448 383-7077	70617-496 383-7648
III 70618-446 383-8737	70618-494 383-9099	70617-449 383-7085	70617-497 383-7655
IV 70618-447 383-8745	70618-495 383-9107	70617-450 383-7093	70617-498 383-7663
V 70618-448 383-8752	70618-496 383-9115	70617-451 383-7101	70617-499 383-7671
VI 70618-449 383-8760	70618-497 383-9123	70617-452 383-7119	70617-500 383-7689

Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Thigh High - Bronze ■ RAL CCL 1 - Closed Toe, Dotted Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-450 383-8778	70618-498 383-9131	70617-453 383-7135	70617-501 383-7697
II 70618-451 383-8786	70618-499 383-9149	70617-454 383-7143	70617-502 383-7705
III 70618-452 383-8794	70618-500 383-9156	70617-455 383-7150	70617-503 383-7713
IV 70618-453 383-8802	70618-501 383-9164	70617-456 383-7168	70617-504 383-7721
V 70618-454 383-8810	70618-502 383-9172	70617-457 383-7176	70617-505 383-7739
VI 70618-455 383-8828	70618-503 383-9180	70617-458 383-7184	70617-506 383-7747

Thigh High - Bronze ■ RAL CCL 1 - Open Toe, Soft Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-420 383-8471	70618-468 383-8836	70617-423 380-7278	70617-471 383-7192
II 70618-421 383-8489	70618-469 383-8844	70617-424 380-7286	70617-472 383-7200
III 70618-422 383-8497	70618-470 383-8851	70617-425 380-7294	70617-473 383-7218
IV 70618-423 383-8505	70618-471 383-8869	70617-426 380-7302	70617-474 383-7226
V 70618-424 383-8513	70618-472 383-8877	70617-427 380-7310	70617-475 383-7234
VI 70618-425 383-8521	70618-473 383-8885	70617-428 380-7328	70617-476 383-7242

Thigh High - Bronze ■ RAL CCL 1 - Closed Toe, Soft Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-426 383-8539	70618-474 383-8893	70617-429 380-7336	70617-477 383-7259
II 70618-427 383-8547	70618-475 383-8901	70617-430 380-7344	70617-478 383-7267
III 70618-428 383-8554	70618-476 383-8919	70617-431 380-7351	70617-479 383-7275
IV 70618-429 383-8562	70618-477 383-8927	70617-432 380-7369	70617-480 383-7283
V 70618-430 383-8570	70618-478 383-8935	70617-433 380-7377	70617-481 383-7291
VI 70618-431 383-8588	70618-479 383-8943	70617-434 380-7385	70617-482 383-7309

Thigh High - Bronze ■ RAL CCL 1 - Open Toe, Lace Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-432 383-8596	70618-480 383-8950	70617-435 380-7393	70617-483 383-7333
II 70618-433 383-8604	70618-481 383-8968	70617-436 380-7401	70617-484 383-7341
III 70618-434 383-8612	70618-482 383-8976	70617-437 380-7419	70617-485 383-7358
IV 70618-435 383-8620	70618-483 383-8984	70617-438 380-7427	70617-486 383-7366
V 70618-436 383-8638	70618-484 383-8992	70617-439 380-7435	70617-487 383-7374
VI 70618-437 383-8646	70618-485 383-9008	70617-440 380-7443	70617-488 383-7382

Thigh High - Bronze ■ RAL CCL 1 - Closed Toe, Lace Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-438 383-8653	70618-486 383-9016	70617-441 380-7450	70617-489 383-7390
II 70618-439 383-8661	70618-487 383-9024	70617-442 380-7468	70617-490 383-7408
III 70618-440 383-8679	70618-488 383-9032	70617-443 380-7476	70617-491 383-7416
IV 70618-441 383-8687	70618-489 383-9040	70617-444 380-7484	70617-492 383-7424
V 70618-442 383-8695	70618-490 383-9057	70617-445 380-7492	70617-493 383-7432
VI 70618-443 383-8703	70618-491 383-9065	70617-446 380-7500	70617-494 383-7440



Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Thigh High - Caramel

RAL CCL 1 - Open Toe, Dotted Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-228 384-2903	70618-360 384-3265	70617-231 384-1186	70617-363 384-1640
II 70618-229 384-2911	70618-361 384-3273	70617-232 384-1194	70617-364 384-1657
III 70618-230 384-2929	70618-362 384-3281	70617-233 384-1202	70617-365 384-1665
IV 70618-231 384-2937	70618-363 384-3299	70617-234 384-1210	70617-366 384-1673
V 70618-232 384-2945	70618-364 384-3307	70617-235 384-1228	70617-367 384-1681
VI 70618-233 384-2952	70618-365 384-3315	70617-236 384-1236	70617-368 384-1699

Thigh High - Caramel

RAL CCL 1 - Closed Toe, Dotted Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-258 384-2960	70618-390 384-3323	70617-261 384-1244	70617-393 384-1731
II 70618-259 384-2978	70618-391 384-3331	70617-262 384-1251	70617-394 384-1749
III 70618-260 384-2986	70618-392 384-3349	70617-263 384-1269	70617-395 384-1756
IV 70618-261 384-2994	70618-393 384-3356	70617-264 384-1277	70617-396 384-1764
V 70618-262 384-3000	70618-394 384-3364	70617-265 384-1285	70617-397 384-1772
VI 70618-263 384-3018	70618-395 384-3372	70617-266 384-1293	70617-398 384-1780

Thigh High - Caramel

RAL CCL 1 - Open Toe, Soft Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-72 384-0949	70618-156 384-3026	70617-72 384-2663	70617-156 384-1335
II 70618-73 384-0956	70618-157 384-3034	70617-73 384-2671	70617-157 384-1343
III 70618-74 384-0964	70618-158 384-3042	70617-74 384-2689	70617-158 384-1350
IV 70618-75 384-0972	70618-159 384-3059	70617-75 384-2697	70617-159 384-1368
V 70618-76 384-0980	70618-160 384-3067	70617-76 384-2705	70617-160 384-1376
VI 70618-77 384-0998	70618-161 384-3075	70617-77 384-2713	70617-161 384-1384

Thigh High - Caramel

RAL CCL 1 - Closed Toe, Soft Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-00 384-2721	70618-84 384-3083	70617-00 384-1004	70617-84 384-1426
II 70618-01 384-2739	70618-85 384-3091	70617-01 384-1012	70617-85 384-1434
III 70618-02 384-2747	70618-86 384-3109	70617-02 384-1020	70617-86 384-1442
IV 70618-03 384-2754	70618-87 384-3117	70617-03 384-1038	70617-87 384-1459
V 70618-04 384-2762	70618-88 384-3125	70617-04 384-1046	70617-88 384-1467
VI 70618-05 384-2770	70618-89 384-3133	70617-05 384-1053	70617-89 384-1475

Thigh High - Caramel

RAL CCL 1 - Open Toe, Lace Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-180 384-2788	70618-198 384-3141	70617-183 384-1061	70617-201 384-1491
II 70618-181 384-2796	70618-199 384-3158	70617-184 384-1079	70617-202 384-1509
III 70618-182 384-2804	70618-200 384-3166	70617-185 384-1087	70617-203 384-1517
IV 70618-183 384-2812	70618-201 384-3174	70617-186 384-1095	70617-204 384-1525
V 70618-184 384-2820	70618-202 384-3182	70617-187 384-1103	70617-205 384-1533
VI 70618-185 384-2838	70618-203 384-3190	70617-188 384-1111	70617-206 384-1541

Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Thigh High - Caramel

RAL CCL 1 - Closed Toe, Lace Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70618-06 384-2846	70618-90 384-3208	70617-06 384-1129	70617-90 384-1574
II 70618-07 384-2853	70618-91 384-3216	70617-07 384-1137	70617-91 384-1582
III 70618-08 384-2861	70618-92 384-3224	70617-08 384-1145	70617-92 384-1590
IV 70618-09 384-2879	70618-93 384-3232	70617-09 384-1152	70617-93 384-1608
V 70618-10 384-2887	70618-94 384-3240	70617-10 384-1160	70617-94 384-1616
VI 70618-11 384-2895	70618-95 384-3257	70617-11 384-1178	70617-95 384-1624

Thigh High - Sand

RAL CCL 2 - Open Toe, Dotted Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70628-126 326-8877 EGD5068	70628-270 326-3142 EGD7947	70627-126 326-2599 EGD4981	70627-270 326-3027 EGD4948
II 70628-127 326-8885 EGD5069	70628-271 326-3159 EGD7948	70627-127 326-2607 EGD4982	70627-271 326-3035 EGD4949
III 70628-128 326-8893 EGD4983	70628-272 326-3167 EGD7949	70627-128 326-2615 EGD4989	70627-272 326-3043 EGD4950
IV 70628-129 326-8901 EGD4978	70628-273 326-3175 EGD7950	70627-129 326-2623 EGD4984	70627-273 326-3050 EGD4951
V 70628-130 326-8919 EGD4979	70628-274 326-3183 EGD7951	70627-130 326-2631 EGD4985	70627-274 326-3068 EGD4958
VI 70628-131 326-8869 EGD4980	70628-275 326-3191 EGD4947	70627-131 326-2649 EGD4986	70627-275 326-3076 EGD4953

Thigh High - Sand

RAL CCL 2 - Closed Toe, Dotted Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70628-150 357-5057 EGD7770	70628-294 357-8150	70627-150 357-5115 EGD7764	70627-294 357-8218
II 70628-151 357-5065 EGD7771	70628-295 357-8168	70627-151 357-5123 EGD7765	70627-295 357-8226
III 70628-152 357-5073 EGD7772	70628-296 357-8176	70627-152 357-5131 EGD7766	70627-296 357-8234
IV 70628-153 357-5081 EGD7773	70628-297 357-8184	70627-153 357-5149 EGD7767	70627-297 357-8242
V 70628-154 357-5099 EGD7774	70628-298 357-8192	70627-154 357-5156 EGD7768	70627-298 357-8259
VI 70628-155 357-5107 EGD7775	70628-299 357-8200	70627-155 357-5164 EGD7769	70627-299 357-8267

Thigh High - Sand

RAL CCL 2 - Open Toe, Soft Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70628-48 358-3499	70628-216 358-3614	70627-48 358-3739	70627-216 358-3861
II 70628-49 358-3507	70628-217 358-3622	70627-49 358-3747	70627-217 358-3879
III 70628-50 358-3515	70628-218 358-3630	70627-50 358-3754	70627-218 358-3887
IV 70628-51 358-3523	70628-219 358-3648	70627-51 358-3762	70627-219 358-3895
V 70628-52 358-3531	70628-220 358-3655	70627-52 358-3770	70627-220 358-3911
VI 70628-53 358-3549	70628-221 358-3663	70627-53 358-3788	70627-221 358-3929

Thigh High - Sand

RAL CCL 2 - Closed Toe, Soft Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 70628-60 357-4431	70628-228 357-4936	70627-60 357-4613	70627-228 357-8036
II 70628-61 357-4449	70628-229 357-4944	70627-61 357-4621	70627-229 357-8044
III 70628-62 357-4456	70628-230 357-4951	70627-62 357-4639	70627-230 357-8051
IV 70628-63 357-4464	70628-231 357-4969	70627-63 357-4647	70627-231 357-8069
V 70628-64 357-4472	70628-232 357-4977	70627-64 357-4654	70627-232 357-8077
VI 70628-65 357-4480	70628-233 357-4985	70627-65 357-4662	70627-233 357-8085



Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Thigh High - Sand ■ RAL CCL 2 - Open Toe, Lace Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-54	357-5875	70628-222	357-7236	70627-54	357-5172	70627-222	357-6626
II	70628-55	357-5883	70628-223	357-7244	70627-55	357-5180	70627-223	357-6634
III	70628-56	357-5891	70628-224	357-7251	70627-56	357-5198	70627-224	357-6642
IV	70628-57	357-5917	70628-225	357-7269	70627-57	357-5206	70627-225	357-6659
V	70628-58	357-5925	70628-226	357-7277	70627-58	357-5214	70627-226	357-6675
VI	70628-59	357-5933	70628-227	357-7285	70627-59	357-5222	70627-227	357-6683

Thigh High - Sand ■ RAL CCL 2 - Closed Toe, Lace Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-66	357-5958	70628-234	357-6733	70627-66	357-6436	70627-234	357-7178
II	70628-67	357-5966	70628-235	357-6741	70627-67	357-6444	70627-235	357-7186
III	70628-68	357-5982	70628-236	357-6766	70627-68	357-6451	70627-236	357-7194
IV	70628-69	357-5990	70628-237	357-6774	70627-69	357-6469	70627-237	357-7202
V	70628-70	357-6006	70628-238	357-6782	70627-70	357-6477	70627-238	357-7210
VI	70628-71	357-6014	70628-239	357-6808	70627-71	357-6485	70627-239	357-7228

Thigh High - Black ■ RAL CCL 2 - Open Toe, Dotted Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-330	326-8786	70628-348	326-3084	70627-330	326-2532	70627-348	326-2656
II	70628-331	326-8794	70628-349	326-3092	70627-331	326-2540	70627-349	326-2664
III	70628-332	326-8810	70628-350	326-3100	70627-332	326-2557	70627-350	326-2672
IV	70628-333	326-8828	70628-351	326-3118	70627-333	326-2565	70627-351	326-2680
V	70628-334	326-8836	70628-352	326-3126	70627-334	326-2573	70627-352	326-2698
VI	70628-335	326-8851	70628-353	326-3134	70627-335	326-2581	70627-353	326-2706

Thigh High - Black ■ RAL CCL 2 - Closed Toe, Dotted Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-156	357-4670	70628-300	357-8093	70627-156	357-4738	70627-300	357-7970
II	70628-157	357-4688	70628-301	357-8101	70627-157	357-4746	70627-301	357-7988
III	70628-158	357-4696	70628-302	357-8119	70627-158	357-4753	70627-302	357-7996
IV	70628-159	357-4704	70628-303	357-8127	70627-159	357-4761	70627-303	357-8002
V	70628-160	357-4712	70628-304	357-8135	70627-160	357-4779	70627-304	357-8010
VI	70628-161	357-4720	70628-305	357-8143	70627-161	357-4787	70627-305	357-8028

Thigh High - Black ■ RAL CCL 2 - Open Toe, Soft Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-312	358-3556	70628-336	358-3671	70627-312	358-3796	70627-336	358-3937
II	70628-313	358-3564	70628-337	358-3689	70627-313	358-3804	70627-337	358-3945
III	70628-314	358-3572	70628-338	358-3697	70627-314	358-3812	70627-338	358-3952
IV	70628-315	358-3580	70628-339	358-3705	70627-315	358-3820	70627-339	358-3960
V	70628-316	358-3598	70628-340	358-3713	70627-316	358-3838	70627-340	358-3978
VI	70628-318	358-3606	70628-341	358-3721	70627-317	358-3846	70627-341	358-3986

Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Thigh High - Black ■ RAL CCL 2 - Closed Toe, Soft Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-72	357-4498	70628-240	357-4993	70627-72	357-4555	70627-240	357-7293
II	70628-73	357-4506	70628-241	357-5008	70627-73	357-4563	70627-241	357-7301
III	70628-74	357-4514	70628-242	357-5016	70627-74	357-4571	70627-242	357-7319
IV	70628-75	357-4522	70628-243	357-5024	70627-75	357-4589	70627-243	357-7327
V	70628-76	357-4530	70628-244	357-5032	70627-76	357-4597	70627-244	357-7335
VI	70628-77	357-4548	70628-245	357-5040	70627-77	357-4605	70627-245	357-7343

Thigh High - Black ■ RAL CCL 2 - Open Toe, Lace Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-317	357-5800	70628-342	357-6493	70627-318	357-5230	70627-342	357-6568
II	70628-319	357-5818	70628-343	357-6501	70627-319	357-5248	70627-343	357-6576
III	70628-320	357-5826	70628-344	357-6527	70627-320	357-5255	70627-344	357-6584
IV	70628-321	357-5842	70628-345	357-6535	70627-321	357-5263	70627-345	357-6592
V	70628-322	357-5859	70628-346	357-6543	70627-322	357-5271	70627-346	357-6600
VI	70628-323	357-5867	70628-347	357-6550	70627-323	357-5289	70627-347	357-6618

Thigh High - Black ■ RAL CCL 2 - Closed Toe, Lace Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-78	357-6022	70628-246	357-6824	70627-78	357-6352	70627-246	357-7350
II	70628-79	357-6030	70628-247	357-6832	70627-79	357-6360	70627-247	357-7368
III	70628-80	357-6048	70628-248	357-6840	70627-80	357-6378	70627-248	357-7376
IV	70628-81	357-6063	70628-249	357-6865	70627-81	357-6394	70627-249	357-7384
V	70628-82	357-6089	70628-250	357-6873	70627-82	357-6402	70627-250	357-7392
VI	70628-83	357-6097	70628-251	357-6881	70627-83	357-6410	70627-251	357-7400

Thigh High - Navy ■ RAL CCL 2 - Open Toe, Dotted Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-378	361-1233	70628-414	361-1290	70627-372	361-1357	70627-414	361-1415
II	70628-379	361-1241	70628-415	361-1308	70627-373	361-1365	70627-415	361-1423
III	70628-380	361-1258	70628-416	361-1316	70627-374	361-1373	70627-416	361-1431
IV	70628-381	361-1266	70628-417	361-1324	70627-375	361-1381	70627-417	361-1449
V	70628-382	361-1274	70628-418	361-1332	70627-376	361-1399	70627-418	361-1456
VI	70628-383	361-1282	70628-419	361-1340	70627-377	361-1407	70627-419	361-1464

Thigh High - Navy ■ RAL CCL 2 - Closed Toe, Dotted Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-144	360-9898	70628-288	360-9971	70627-144	361-0037	70627-288	361-0094
II	70628-145	360-9914	70628-289	360-9989	70627-145	361-0045	70627-289	361-0102
III	70628-146	360-9922	70628-290	360-9997	70627-146	361-0052	70627-290	361-0110
IV	70628-147	360-9930	70628-291	361-0003	70627-147	361-0060	70627-291	361-0128
V	70628-148	360-9955	70628-292	361-0011	70627-148	361-0078	70627-292	361-0136
VI	70628-149	360-9963	70628-293	361-0029	70627-149	361-0086	70627-293	361-0144



Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Thigh High - Navy ■ RAL CCL 2 - Open Toe, Soft Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-366	361-0151	70628-402	361-0219	70627-360	361-0276	70627-402	361-0334
II	70628-367	361-0169	70628-403	361-0227	70627-361	361-0284	70627-403	361-0342
III	70628-368	361-0177	70628-404	361-0235	70627-362	361-0292	70627-404	361-0359
IV	70628-369	361-0185	70628-405	361-0243	70627-363	361-0300	70627-405	361-0367
V	70628-370	361-0193	70628-406	361-0250	70627-364	361-0318	70627-406	361-0375
VI	70628-371	361-0201	70628-407	361-0268	70627-365	361-0326	70627-407	361-0383

Thigh High - Navy ■ RAL CCL 2 - Closed Toe, Soft Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-36	360-9344	70628-204	360-9476	70627-36	360-9534	70627-204	360-9591
II	70628-37	360-9351	70628-205	360-9484	70627-37	360-9542	70627-205	360-9609
III	70628-38	360-9369	70628-206	360-9492	70627-38	360-9559	70627-206	360-9617
IV	70628-39	360-9377	70628-207	360-9500	70627-39	360-9567	70627-207	360-9625
V	70628-40	360-9385	70628-208	360-9518	70627-40	360-9575	70627-208	360-9633
VI	70628-41	360-9393	70628-209	360-9526	70627-41	360-9583	70627-209	360-9641

Thigh High - Navy ■ RAL CCL 2 - Open Toe, Lace Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-372	361-0391	70628-408	361-0458	70627-366	361-0516	70627-408	361-0573
II	70628-373	361-0409	70628-409	361-0466	70627-367	361-0524	70627-409	361-0581
III	70628-374	361-0417	70628-410	361-0474	70627-368	361-0532	70627-410	361-0599
IV	70628-375	361-0425	70628-411	361-0482	70627-369	361-0540	70627-411	361-0607
V	70628-376	361-0433	70628-412	361-0490	70627-370	361-0557	70627-412	361-0615
VI	70628-377	361-0441	70628-413	361-0508	70627-371	361-0565	70627-413	361-0623

Thigh High - Navy ■ RAL CCL 2 - Closed Toe, Lace Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-42	360-9658	70628-210	360-9716	70627-42	360-9773	70627-210	360-9831
II	70628-43	360-9666	70628-211	360-9724	70627-43	360-9781	70627-211	360-9849
III	70628-44	360-9674	70628-212	360-9732	70627-44	360-9799	70627-212	360-9856
IV	70628-45	360-9682	70628-213	360-9740	70627-45	360-9807	70627-213	360-9864
V	70628-46	360-9690	70628-214	360-9757	70627-46	360-9815	70627-214	360-9872
VI	70628-47	360-9708	70628-215	360-9765	70627-47	360-9823	70627-215	360-9880

Thigh High - Bronze ■ RAL CCL 2 - Open Toe, Dotted Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-444	383-9446	70628-492	383-9818	70627-444	383-7994	70627-492	383-8356
II	70628-445	383-9453	70628-493	383-9826	70627-445	383-8000	70627-493	383-8364
III	70628-446	383-9461	70628-494	383-9834	70627-446	383-8018	70627-494	383-8372
IV	70628-447	383-9479	70628-495	383-9842	70627-447	383-8026	70627-495	383-8380
V	70628-448	383-9487	70628-496	383-9859	70627-448	383-8034	70627-496	383-8398
VI	70628-449	383-9495	70628-497	383-9867	70627-449	383-8042	70627-497	383-8406

Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Thigh High - Bronze ■ RAL CCL 2 - Closed Toe, Dotted Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-450	383-9511	70628-498	383-9875	70627-450	383-8059	70627-498	383-8414
II	70628-451	383-9529	70628-499	383-9883	70627-451	383-8067	70627-499	383-8422
III	70628-452	383-9537	70628-500	383-9891	70627-452	383-8075	70627-500	383-8430
IV	70628-453	383-9545	70628-501	383-9909	70627-453	383-8083	70627-501	383-8448
V	70628-454	383-9552	70628-502	383-9917	70627-454	383-8091	70627-502	383-8455
VI	70628-455	383-9560	70628-503	383-9925	70627-455	383-8109	70627-503	383-8463

Thigh High - Bronze ■ RAL CCL 2 - Open Toe, Soft Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-420	383-9198	70628-468	383-9578	70627-420	383-7754	70627-468	383-8117
II	70628-421	383-9198	70628-469	383-9586	70627-421	383-7762	70627-469	383-8125
III	70628-422	383-9214	70628-470	383-9594	70627-422	383-7770	70627-470	383-8133
IV	70628-423	383-9222	70628-471	383-9602	70627-423	383-7788	70627-471	383-8141
V	70628-424	383-9230	70628-472	383-9610	70627-424	383-7796	70627-472	383-8158
VI	70628-425	383-9248	70628-473	383-9628	70627-425	383-7804	70627-473	383-8166

Thigh High - Bronze ■ RAL CCL 2 - Closed Toe, Soft Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-426	383-9263	70628-474	383-9636	70627-426	383-7812	70627-474	383-8174
II	70628-427	383-9271	70628-475	383-9644	70627-427	383-7820	70627-475	383-8182
III	70628-428	383-9289	70628-476	383-9651	70627-428	383-7838	70627-476	383-8190
IV	70628-429	383-9297	70628-477	383-9669	70627-429	383-7846	70627-477	383-8208
V	70628-430	383-9305	70628-478	383-9677	70627-430	383-7853	70627-478	383-8216
VI	70628-431	383-9313	70628-479	383-9685	70627-431	383-7861	70627-479	383-8224

Thigh High - Bronze ■ RAL CCL 2 - Open Toe, Lace Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-432	383-9321	70628-480	383-9693	70627-432	383-7879	70627-480	383-8232
II	70628-433	383-9339	70628-481	383-9701	70627-433	383-7887	70627-481	383-8240
III	70628-434	383-9347	70628-482	383-9719	70627-434	383-7895	70627-482	383-8257
IV	70628-435	383-9354	70628-483	383-9727	70627-435	383-7903	70627-483	383-8265
V	70628-436	383-9362	70628-484	383-9735	70627-436	383-7911	70627-484	383-8273
VI	70628-437	383-9370	70628-485	383-9743	70627-437	383-7929	70627-485	383-8281

Thigh High - Bronze ■ RAL CCL 2 - Closed Toe, Lace Silicone Band

	Petite		Petite, wider circumference		Regular		Regular, wider circumference	
I	70628-438	383-9388	70628-486	383-9750	70627-438	383-7937	70627-486	383-8299
II	70628-439	383-9396	70628-487	383-9768	70627-439	383-7945	70627-487	383-8307
III	70628-440	383-9404	70628-488	383-9776	70627-440	383-7952	70627-488	383-8315
IV	70628-441	383-9412	70628-489	383-9784	70627-441	383-7960	70627-489	383-8323
V	70628-442	383-9420	70628-490	383-9792	70627-442	383-7978	70627-490	383-8331
VI	70628-443	383-9438	70628-491	383-9800	70627-443	383-7986	70627-491	383-8349



Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Tights - Caramel

RAL CCL 1 - Open Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

RAL CCL 1 - Closed Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

Tights - Sand

RAL CCL 2 - Open Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

RAL CCL 2 - Closed Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

Tights - Black

RAL CCL 2 - Open Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

RAL CCL 2 - Closed Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

Tights - Navy

RAL CCL 2 - Open Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

RAL CCL 2 - Closed Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

Tights - Bronze

RAL CCL 2 - Open Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

RAL CCL 2 - Closed Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Tights - Caramel

RAL CCL 2 - Open Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

RAL CCL 2 - Closed Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

Tights for Men - Sand

RAL CCL 1 - Open Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

RAL CCL 1 - Closed Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

Tights for Men - Sand

RAL CCL 2 - Open Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

RAL CCL 2 - Closed Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

Tights for Men - Black

RAL CCL 1 - Open Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

RAL CCL 1 - Closed Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

Tights for Men - Black

RAL CCL 2 - Open Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.

RAL CCL 2 - Closed Toe

Table with columns: Petite, Regular. Rows I-VI for BSN codes and sizes.



Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Chap Style - Sand ■

RAL CCL 1 - Open Toe

Left, petite		Left, regular		Right, petite		Right, regular	
I	70620-54 EGD7692	70619-54	EGD7698	70620-24	EGD7704	70619-24	EGD7710
II	70620-55 EGD7693	70619-55	EGD7699	70620-25	EGD7705	70619-25	EGD7711
III	70620-56 EGD7694	70619-56	EGD7700	70620-26	EGD7706	70619-26	EGD7712
IV	70620-57 EGD7695	70619-57	EGD7701	70620-27	EGD7707	70619-27	EGD7713
V	70620-58 EGD7696	70619-58	EGD7702	70620-28	EGD7708	70619-28	EGD7714
VI	70620-59 EGD7697	70619-59	EGD7703	70620-29	EGD7709	70619-29	EGD7715

Chap Style - Sand ■

RAL CCL 1 - Closed Toe

Left, petite		Left, regular		Right, petite		Right, regular	
I	70620-144	70619-144		70620-114	70619-114		
II	70620-145	70619-145		70620-115	70619-115		
III	70620-146	70619-146		70620-116	70619-116		
IV	70620-147	70619-147		70620-117	70619-117		
V	70620-148	70619-148		70620-118	70619-118		
VI	70620-149	70619-149		70620-119	70619-119		

Chap Style - Sand ■

RAL CCL 1 - Open Toe

Pair, petite		Pair, regular		Pair, petite		Pair, regular	
I	70620-84	70619-84		70620-174	70619-174		
II	70620-85	70619-85		70620-175	70619-175		
III	70620-86	70619-86		70620-176	70619-176		
IV	70620-87	70619-87		70620-177	70619-177		
V	70620-88	70619-88		70620-178	70619-178		
VI	70620-89	70619-89		70620-179	70619-179		

RAL CCL 1 - Closed Toe

Chap Style - Sand ■

RAL CCL 2 - Open Toe

Left, petite		Left, regular		Right, petite		Right, regular	
I	70630-42 EGD7782	70629-42	EGD7776	70630-36	EGD7794	70629-36	EGD7788
II	70630-43 EGD7783	70629-43	EGD7777	70630-37	EGD7795	70629-37	EGD7789
III	70630-44 EGD7784	70629-44	EGD7778	70630-38	EGD7796	70629-38	EGD7790
IV	70630-45 EGD7785	70629-45	EGD7779	70630-39	EGD7797	70629-39	EGD7791
V	70630-46 EGD7786	70629-46	EGD7780	70630-40	EGD7798	70629-40	EGD7792
VI	70630-47 EGD7787	70629-47	EGD7781	70630-41	EGD7799	70629-41	EGD7793

Chap Style - Sand ■

RAL CCL 2 - Closed Toe

Left, petite		Left, regular		Right, petite		Right, regular	
I	70630-60	70629-60		70630-54	70629-54		
II	70630-61	70629-61		70630-55	70629-55		
III	70630-62	70629-62		70630-56	70629-56		
IV	70630-63	70629-63		70630-57	70629-57		
V	70630-64	70629-64		70630-58	70629-58		
VI	70630-65	70629-65		70630-59	70629-59		

Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Chap Style - Sand ■

RAL CCL 2 - Open Toe

Pair, petite		Pair, regular		Pair, petite		Pair, regular	
I	70630-48	70629-48		70630-66	70629-66		
II	70630-49	70629-49		70630-67	70629-67		
III	70630-50	70629-50		70630-68	70629-68		
IV	70630-51	70629-51		70630-69	70629-69		
V	70630-52	70629-52		70630-70	70629-70		
VI	70630-53	70629-53		70630-71	70629-71		

RAL CCL 2 - Closed Toe

Maternity - Sand ■

RAL CCL 1 - Open Toe

Petite		Regular		Petite		Regular	
I	70622-102	70621-102		70622-42	70621-42		
II	70622-103	70621-103		70622-43	70621-43		
III	70622-104	70621-104		70622-44	70621-44		
IV	70622-105	70621-105		70622-45	70621-45		
V	70622-106	70621-106		70622-46	70621-46		
VI	70622-107	70621-107		70622-47	70621-47		

RAL CCL 1 - Closed Toe

Maternity - Black ■

RAL CCL 1 - Open Toe

Petite		Regular		Petite		Regular	
I	70622-96	70621-96		70622-54	70621-54		
II	70622-97	70621-97		70622-55	70621-55		
III	70622-98	70621-98		70622-56	70621-56		
IV	70622-99	70621-99		70622-57	70621-57		
V	70622-100	70621-100		70622-58	70621-58		
VI	70622-101	70621-101		70622-59	70621-59		

RAL CCL 1 - Closed Toe

Maternity - Sand ■

RAL CCL 2 - Open Toe

Petite		Regular		Petite		Regular	
I	70632-60	70631-60		70632-66	70631-66		
II	70632-61	70631-61		70632-67	70631-67		
III	70632-62	70631-62		70632-68	70631-68		
IV	70632-63	70631-63		70632-69	70631-69		
V	70632-64	70631-64		70632-70	70631-70		
VI	70632-65	70631-65		70632-71	70631-71		

RAL CCL 2 - Closed Toe

Maternity - Black ■

RAL CCL 2 - Open Toe

Petite		Regular		Petite		Regular	
I	70632-120	70631-120		70632-78	70631-78		
II	70632-121	70631-121		70632-79	70631-79		
III	70632-122	70631-122		70632-80	70631-80		
IV	70632-123	70631-123		70632-81	70631-81		
V	70632-124	70631-124		70632-82	70631-82		
VI	70632-125	70631-125		70632-83	70631-83		

RAL CCL 2 - Closed Toe



Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Knee High - Anthracite ■ RAL CCL 1 - Closed Toe

	Petite		Regular	
I	73576-30	389-1546	73575-30	389-1538
II	73576-31	389-2049	73575-31	389-1694
III	73576-32	389-2056	73575-32	389-1702
IV	73576-33	389-2064	73575-33	389-1710
V	73576-34	389-2072	73575-34	389-1728
VI	73576-35	389-2080	73575-35	389-1736

Knee High - Black ■ RAL CCL 1 - Closed Toe

	Petite		Regular	
	73576-24	389-2098	73575-24	389-1744 EGD4973
	73576-25	389-2106	73575-25	389-1751 EGD4974
	73576-26	389-2114	73575-26	389-1769 EGD4891
	73576-27	389-2122	73575-27	389-1777 EGD4936
	73576-28	389-2130	73575-28	389-1785 EGD4937
	73576-29	389-2148	73575-29	389-1793 EGD4938

Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Knee High - Black ■ RAL CCL 2 - Open Toe

	Petite		Regular	
I	73584-52	389-4763	73583-48	389-4581
II	73584-48	389-4771	73583-49	389-4599
III	73584-49	389-4789	73583-50	389-4607
IV	73584-50	389-4797	73583-51	389-4615
V	73584-51	389-4805	73583-52	389-4623
VI	73584-53	389-4813	73583-53	389-4631

RAL CCL 2 - Closed Toe

	Petite		Regular	
	73584-36	389-3633	73583-36	389-3278
	73584-37	389-3641	73583-37	389-3286
	73584-38	389-3658	73583-38	389-3294
	73584-39	389-3666	73583-39	389-3302
	73584-40	389-3674	73583-40	389-3310
	73584-41	389-3682	73583-41	389-3328

Knee High - Bronze ■ RAL CCL 1 - Closed Toe

	Petite		Regular	
I	73576-96	389-2155	73575-84	389-1801
II	73576-97	389-2163	73575-85	389-1819
III	73576-98	389-2171	73575-86	389-1827
IV	73576-99	389-2189	73575-87	389-1835
V	73576-100	389-2197	73575-88	389-1843
VI	73576-101	389-2205	73575-89	389-1850

Knee High - Café ■ RAL CCL 1 - Closed Toe

	Petite		Regular	
	73576-12	389-2213	73575-12	389-1868
	73576-13	389-2221	73575-13	389-1876
	73576-14	389-2239	73575-14	389-1884
	73576-15	389-2247	73575-15	389-1892
	73576-16	389-2254	73575-16	389-1900
	73576-17	389-2262	73575-17	389-1918

Knee High - Café ■ RAL CCL 2 - Closed Toe

	Petite		Regular	
I	73584-18	389-3757	73583-18	389-3393
II	73584-19	389-3765	73583-19	389-3401
III	73584-20	389-3773	73583-20	389-3419
IV	73584-21	389-3781	73583-21	389-3427
V	73584-22	389-3799	73583-22	389-3435
VI	73584-23	389-3807	73583-23	389-3443

Knee High - Caramel ■ RAL CCL 1 - Closed Toe

	Petite		Regular	
I	73576-00	389-2270	73575-00	389-1926
II	73576-01	389-2288	73575-01	389-1934
III	73576-02	389-2296	73575-02	389-1942
IV	73576-03	389-2304	73575-03	389-1959
V	73576-04	389-2312	73575-04	389-1967
VI	73576-05	389-2320	73575-05	389-1975

Knee High - Natural ■ RAL CCL 1 - Closed Toe

	Petite		Regular	
	73576-18	389-2338	73575-18	389-1983 EGD4968
	73576-19	389-2346	73575-19	389-1991 EGD4969
	73576-20	389-2353	73575-20	389-2007 EGD4975
	73576-21	389-2361	73575-21	389-2015 EGD4976
	73576-22	389-2379	73575-22	389-2023 EGD4971
	73576-23	389-2387	73575-23	389-2031 EGD4972

Knee High - Caramel ■ RAL CCL 2 - Open Toe

	Petite		Regular	
I	73584-00	389-4821	73583-00	389-4649
II	73584-01	389-4839	73583-01	389-4656
III	73584-02	389-4847	73583-02	389-4664
IV	73584-03	389-4854	73583-03	389-4672
V	73584-04	389-4862	73583-04	389-4680
VI	73584-05	389-4870	73583-05	389-4698

RAL CCL 2 - Closed Toe

	Petite		Regular	
	73584-06	389-3815	73583-06	389-3450
	73584-07	389-3823	73583-07	389-3468
	73584-08	389-3831	73583-08	389-3476
	73584-09	389-3849	73583-09	389-3484
	73584-10	389-3856	73583-10	389-3492
	73584-11	389-3864	73583-11	389-3500

Knee High - Anthracite ■ RAL CCL 2 - Closed Toe

	Petite		Regular	
I	73584-42	389-3575	73583-42	389-3211
II	73584-43	389-3583	73583-43	389-3229
III	73584-44	389-3591	73583-44	389-3237
IV	73584-45	389-3609	73583-45	389-3245
V	73584-46	389-3617	73583-46	389-3252
VI	73584-47	389-3625	73583-47	389-3260

Knee High - Bronze ■ RAL CCL 2 - Closed Toe

	Petite		Regular	
	73584-115	389-3690	73583-102	389-3336
	73584-116	389-3708	73583-103	389-3344
	73584-117	389-3716	73583-104	389-3351
	73584-118	389-3724	73583-105	389-3369
	73584-119	389-3732	73583-106	389-3377
	73584-120	389-3740	73583-107	389-3385

Knee High - Natural ■ RAL CCL 2 - Open Toe

	Petite		Regular	
I	73584-24	389-4888	73583-24	389-4706 EGD4918
II	73584-25	389-4896	73583-25	389-4714 EGD4919
III	73584-26	389-4904	73583-26	389-4722 EGD4939
IV	73584-27	389-4912	73583-27	389-4730 EGD4940
V	73584-28	389-4920	73583-28	389-4748 EGD4920
VI	73584-29	389-4938	73583-29	389-4755 EGD4977

RAL CCL 2 - Closed Toe

	Petite		Regular	
	73584-30	389-3872	73583-30	389-3518
	73584-31	389-3880	73583-31	389-3526
	73584-32	389-3898	73583-32	389-3534
	73584-33	389-3906	73583-33	389-3542
	73584-34	389-3914	73583-34	389-3559
	73584-35	389-3922	73583-35	389-3567



Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Thigh High - Caramel

RAL CCL 2 - Open Toe, Decorative Dotted Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 73586-210 390-8076	73586-258 390-8258	73585-202 390-6633	73585-250 390-7177
II 73586-211 390-8084	73586-259 390-8266	73585-203 390-6641	73585-251 390-7185
III 73586-212 390-8092	73586-260 390-8274	73585-204 390-6658	73585-252 390-7193
IV 73586-213 390-8100	73586-261 390-8282	73585-205 390-6666	73585-253 390-7201
V 73586-214 390-8118	73586-262 390-8290	73585-206 390-6674	73585-254 390-7219
VI 73586-215 390-8126	73586-263 390-8308	73585-207 390-6682	73585-255 390-7227

Thigh High - Caramel

RAL CCL 2 - Closed Toe, Decorative Dotted Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 73586-222 390-7532	73586-270 390-7896	73585-214 390-6450	73585-262 390-6997
II 73586-223 390-7540	73586-271 390-7904	73585-215 390-6468	73585-263 390-7003
III 73586-224 390-7557	73586-272 390-7912	73585-216 390-6476	73585-264 390-7011
IV 73586-225 390-7565	73586-273 390-7920	73585-217 390-6484	73585-265 390-7029
V 73586-226 390-7573	73586-274 390-7938	73585-218 390-6492	73585-266 390-7037
VI 73586-227 390-7581	73586-275 390-7946	73585-219 390-6500	73585-267 390-7045

Thigh High - Natural

RAL CCL 2 - Open Toe, Decorative Dotted Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 73586-216 390-8134	73586-264 390-8316	73585-208 390-6690	73585-256 390-7235
II 73586-217 390-8142	73586-265 390-8324	73585-209 390-6708	73585-257 390-7243
III 73586-218 390-8159	73586-266 390-8332	73585-210 390-6716	73585-258 390-7250
IV 73586-219 390-8167	73586-267 390-8340	73585-211 390-6724	73585-259 390-7268
V 73586-220 390-8175	73586-268 390-8357	73585-212 390-6732	73585-260 390-7276
VI 73586-221 390-8183	73586-269 390-8365	73585-213 390-6740	73585-261 390-7284

Thigh High - Natural

RAL CCL 2 - Closed Toe, Decorative Dotted Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 73586-240 390-7599	73586-288 390-7953	73585-232 390-6518	73585-280 390-7052
II 73586-241 390-7607	73586-289 390-7961	73585-233 390-6526	73585-281 390-7060
III 73586-242 390-7615	73586-290 390-7979	73585-234 390-6534	73585-282 390-7078
IV 73586-243 390-7623	73586-291 390-7987	73585-235 390-6542	73585-283 390-7086
V 73586-244 390-7631	73586-292 390-7995	73585-236 390-6559	73585-284 390-7094
VI 73586-245 390-7649	73586-293 390-8001	73585-237 390-6567	73585-285 390-7102

Thigh High - Anthracite

RAL CCL 2 - Closed Toe, Soft Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 73586-84 390-4984	73586-198 390-5494	73585-84 390-3366	73585-190 390-3721
II 73586-85 390-4992	73586-199 390-5502	73585-85 390-3374	73585-191 390-3739
III 73586-86 390-5007	73586-200 390-5510	73585-86 390-3382	73585-192 390-3747
IV 73586-87 390-5015	73586-201 390-5528	73585-87 390-3390	73585-193 390-3754
V 73586-88 390-5023	73586-202 390-5536	73585-88 390-3408	73585-194 390-3762
VI 73586-89 390-5031	73586-203 390-5544	73585-89 390-3416	73585-195 390-3770

Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Thigh High - Black

RAL CCL 2 - Open Toe, Soft Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 73586-306 390-5858	73586-314 390-6039	73585-298 390-4083	73585-306 390-4273
II 73586-110 390-5866	73586-315 390-6047	73585-102 390-4091	73585-307 390-4281
III 73586-111 390-5874	73586-316 390-6054	73585-103 390-4109	73585-308 390-4299
IV 73586-112 390-5882	73586-317 390-6054	73585-104 390-4117	73585-309 390-4307
V 73586-113 390-5890	73586-318 390-6070	73585-105 390-4125	73585-310 390-4315
VI 73586-307 390-5908	73586-319 390-6088	73585-299 390-4133	73585-311 390-4323

Thigh High - Black

RAL CCL 2 - Closed Toe, Soft Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 73586-72 390-5049	73586-186 390-5551	73585-72 390-3424	73585-178 390-3788
II 73586-73 390-5056	73586-187 390-5569	73585-73 390-3432	73585-179 390-3796
III 73586-74 390-5064	73586-188 390-5577	73585-74 390-3440	73585-180 390-3804
IV 73586-75 390-5072	73586-189 390-5585	73585-75 390-3457	73585-181 390-3812
V 73586-76 390-5080	73586-190 390-5593	73585-76 390-3465	73585-182 390-3820
VI 73586-77 390-5098	73586-191 390-5601	73585-77 390-3473	73585-183 390-3838

Thigh High - Bronze

RAL CCL 2 - Closed Toe, Soft Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 73586-457 390-5106	73586-469 390-5619	73585-449 390-3481	73585-461 390-3846
II 73586-458 390-5114	73586-470 390-5627	73585-450 390-3499	73585-462 390-3853
III 73586-459 390-5122	73586-471 390-5635	73585-451 390-3507	73585-463 390-3861
IV 73586-460 390-5130	73586-472 390-5643	73585-452 390-3515	73585-464 390-3879
V 73586-461 390-5148	73586-473 390-5650	73585-453 390-3523	73585-465 390-3887
VI 73586-462 390-5155	73586-474 390-5668	73585-454 390-3531	73585-466 390-3895

Thigh High - Café

RAL CCL 2 - Closed Toe, Soft Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 73586-36 390-5163	73586-150 390-5676	73585-36 390-3549	73585-142 390-3903
II 73586-37 390-5171	73586-151 390-5684	73585-37 390-3556	73585-143 390-3911
III 73586-38 390-5189	73586-152 390-5692	73585-38 390-3564	73585-144 390-3929
IV 73586-39 390-5197	73586-153 390-5700	73585-39 390-3572	73585-145 390-3937
V 73586-40 390-5205	73586-154 390-5718	73585-40 390-3580	73585-146 390-3945
VI 73586-41 390-5213	73586-155 390-5726	73585-41 390-3598	73585-147 390-3952

Thigh High - Caramel

RAL CCL 2 - Open Toe, Soft Silicone Band

Petite	Petite, wider circumference	Regular	Regular, wider circumference
I 73586-00 390-5916	73586-114 390-6096	73585-00 390-4141	73585-106 390-4331
II 73586-01 390-5924	73586-115 390-6104	73585-01 390-4158	73585-107 390-4349
III 73586-02 390-5932	73586-116 390-6112	73585-02 390-4166	73585-108 390-4356
IV 73586-03 390-5940	73586-117 390-6120	73585-03 390-4174	73585-109 390-4364
V 73586-04 390-5957	73586-118 390-6138	73585-04 390-4182	73585-110 390-4372
VI 73586-05 390-5965	73586-119 390-6146	73585-05 390-4190	73585-111 390-4380



Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Maternity Tights - Anthracite ■ RAL CCL 1 - Closed Toe

	Petite	Regular
I	73582-66	73581-66
II	73582-67	73581-67
III	73582-68	73581-68
IV	73582-69	73581-69
V	73582-70	73581-70
VI	73582-71	73581-71

Maternity Tights - Black ■ RAL CCL 1 - Closed Toe

	Petite	Regular
I	73582-54	73581-54
II	73582-55	73581-55
III	73582-56	73581-56
IV	73582-57	73581-57
V	73582-58	73581-58
VI	73582-59	73581-59

Maternity Tights - Bronze ■ RAL CCL 1 - Closed Toe

	Petite	Regular
I	73582-222	73581-210
II	73582-223	73581-211
III	73582-224	73581-212
IV	73582-225	73581-213
V	73582-226	73581-214
VI	73582-227	73581-215

Maternity Tights - Café ■ RAL CCL 1 - Closed Toe

	Petite	Regular
I	73582-30	73581-30
II	73582-31	73581-31
III	73582-32	73581-32
IV	73582-33	73581-33
V	73582-34	73581-34
VI	73582-35	73581-35

Maternity Tights - Caramel ■ RAL CCL 1 - Closed Toe

	Petite	Regular
I	73582-06	73581-06
II	73582-07	73581-07
III	73582-08	73581-08
IV	73582-09	73581-09
V	73582-10	73581-10
VI	73582-11	73581-11

Maternity Tights - Natural ■ RAL CCL 1 - Closed Toe

	Petite	Regular
I	73582-42	73581-42
II	73582-43	73581-43
III	73582-44	73581-44
IV	73582-45	73581-45
V	73582-46	73581-46
VI	73582-47	73581-47

Maternity Tights - Anthracite ■ RAL CCL 2 - Closed Toe

	Petite	Regular
I	73590-90	73589-90
II	73590-91	73589-91
III	73590-92	73589-92
IV	73590-93	73589-93
V	73590-94	73589-94
VI	73590-95	73589-95

Maternity Tights - Bronze ■ RAL CCL 2 - Closed Toe

	Petite	Regular
I	73590-324	73589-312
II	73590-325	73589-313
III	73590-326	73589-314
IV	73590-327	73589-315
V	73590-328	73589-316
VI	73590-329	73589-317

Maternity Tights - Black ■ RAL CCL 2 - Open Toe

	Petite	Regular
I	73590-102	73589-102
II	73590-103	73589-103
III	73590-104	73589-104
IV	73590-105	73589-105
V	73590-106	73589-106
VI	73590-107	73589-107

Maternity Tights - Black ■ RAL CCL 2 - Closed Toe

	Petite	Regular
I	73590-78	73589-78
II	73590-79	73589-79
III	73590-80	73589-80
IV	73590-81	73589-81
V	73590-82	73589-82
VI	73590-83	73589-83

Ordering Information

Size BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC BSN Code PIP NHSSC

Maternity Tights - Café ■ RAL CCL 2 - Closed Toe

	Petite	Regular
I	73590-42	73589-42
II	73590-43	73589-43
III	73590-44	73589-44
IV	73590-45	73589-45
V	73590-46	73589-46
VI	73590-47	73589-47

Maternity Tights - Caramel ■ RAL CCL 2 - Open Toe

	Petite	Regular
I	73590-12	73589-12
II	73590-13	73589-13
III	73590-14	73589-14
IV	73590-15	73589-15
V	73590-16	73589-16
VI	73590-17	73589-17

Maternity Tights - Caramel ■ RAL CCL 2 - Closed Toe

	Petite	Regular
I	73590-18	73589-18
II	73590-19	73589-19
III	73590-20	73589-20
IV	73590-21	73589-21
V	73590-22	73589-22
VI	73590-23	73589-23

Maternity Tights - Natural ■ RAL CCL 2 - Open Toe

	Petite	Regular
I	73590-60	73589-60
II	73590-61	73589-61
III	73590-62	73589-62
IV	73590-63	73589-63
V	73590-64	73589-64
VI	73590-65	73589-65


Maternity Tights - Natural ■ RAL CCL 2 - Closed Toe

	Petite	Regular
I	73590-66	73589-66
II	73590-67	73589-67
III	73590-68	73589-68
IV	73590-69	73589-69
V	73590-70	73589-70
VI	73590-71	73589-71

JOBST® UltraSheer and JOBST® Opaque Measuring Guide

RAL Compression classes 1 (18-21mmHg) and 2 (23-32mmHg)						
Measurement Point	I	II	III	IV	V	VI
cG silicone band	43-48	45-52	49-56	53-60	56-64	60-68
cG silicone band wide	49-57	53-62	57-67	61-72	65-77	69-81
cG	43-57	45-62	49-67	53-72	56-77	60-81
cF	39-52	41-56	44-60	47-65	50-69	53-73
cE	30-37	33-40	35-43	37-45	39-48	41-51
cD	27-33	29-36	32-39	34-42	36-45	38-48
cC	28-34	30-37	33-40	35-43	37-46	39-49
cB1	23-27	24-29	26-32	29-35	31-37	33-39
cB	18-20	20-22	22-24	24-26	26-28	28-30
cY	26-31	28-33	29-35	31-37	32-38	33-40
cA	17-22	19-24	21-26	23-29	25-32	27-34

Length measurements for all sizes (in cm):
AD: ID regular 39-45, ID petite 34-39 AG: IG regular 70-80, IG petite 60-69 AT/AG-T: IK regular 75-85, IK petite 65-75



Information



Customer Service

Our experienced UK Customer Service team use their expertise to process orders efficiently and answer product queries.

- Call 0845 122 3600
- Direct fax 0845 122 3450
- Email: compression.uk@bsnmedical.com



Training and Education

Through the JOBST Academy you can create your own, local, bespoke study day. We have a range of therapy focused training modules covering both theory and practical sessions.

- Contact your local BSN medical Compression Therapy Account Manager for more details



Clinical Support

Our free compression therapy helpline is manned by our experienced Clinical Nurse Specialist who is available to share her expertise with you. She will answer your clinical queries and offer treatment advice from bandaging through to measuring for compression garments.

- Call 0800 389 8424



Compression Therapy Account Manager

Your local BSN medical Compression Therapy Account Manager is your first point of contact for our entire range of compression therapy products and as such is able to offer you advice, support and detailed product information.

- Call your local BSN medical Compression Therapy Account Manager directly or call 0845 122 3600 to get further contact details