

USCOPE SPECIFICATIONS

UE3011


Optical system	
Field of view	120°
Line of vision	0°
Visible range	3-50mm
Illumination	Optical fiber
Insertion portion	
Bending section1	>270° up, >270° down
Diameter of insertion cord	≤ 9F
Maximum diameter of distal tip	9F
Maximum diameter of insertion portion	9F
Working length of insertion cord	630mm
Channel	
Average inner diameter	3.6F
Minimum inner diameter of working channel	3.6F

IMAGE SYSTEM SPECIFICATIONS

UTV100

Display	
Max. resolution	1,024 * 768
Orientation	Landscape
Display type	12.1" color TFT LCD
Startup time	Approximately 20 seconds
Storage capacity	8 GB
Dimensions	
Width	315mm
Height	308mm
Thickness	187mm
Weight	4.3 kg
Connections	
USB connection	Type A
AV connection	RCA connector
HDMI connection	Optional

Zhuhai Pusen Medical Technology Co., Ltd.
Address: 5/F Building 1, No 33, Keji San Road, High-tech Zone,
Tangjiaowan Town, Zhuhai, Guangdong, PRC.
Tel: 0086 756 6880096 Fax: 0086 756 6880401
<http://www.pusenmedical.com> Email: Pusen@pusenmedical.com



Note: Specifications indicated in the above table are subject to change without any notice.



SINGLE-USE URETEROSCOPE UE3011

WE BRING IDEAL FOR LITHOTRIPSY

NO-CROSS INFECTION
NO-MAINTENANCES
NO-STERILIZATION
ALL IN ONE DESIGN



PUSEN

PUSEN

HIGH IMAGE QUALITY

No Distortion when laser is On



MONITOR



FIBER


VS




UE3011

IMPROVED HANDLING

LOCATION LOCK




FLEXIBLE DESIGN



With 630mm Working length

WIDE ANGULATION



SLIM DESIGN



WHEEL CONTROL



HANDLE LOCATION LEFT OR RIGHT

