



### Oscimed Ergo II Plaster Saw

#### Ergonomic, compact lightweight and easy to use

The Oscimed ERGO II is a compact, lightweight cast saw ergonomically shaped to enable precision work without any effort. The ERGO II is robustly built and equipped with a cooling fan for intensive use without overheating. The device is powered by a safe 24 volts and complies with all the latest medical standards.

**Power:** 130 Watt – 24VDC  
**Length:** 230mm  
**Diameter:** 58mm  
**Weight:** 860g  
**Blade Oscillation:** 14,500  
**Noise Level:** 67 dBA  
**Standard:** EC conformity EN 60601  
**Other:** Instant stop button



### Oscimed 2000 Plaster Saw

#### The Oscimed 2000 is a lightweight, low noise ergonomic saw

It is extensively used throughout the world and is recommended for its ease of use. With a 24 volt brushless motor it is very safe to use and incorporates instant blade change without the use of tools.



#### Key Features

- Saws are ergonomically designed to be light and can be used single handed.
- Blades are designed to be changed instantly.
- Blades have a noise suppressing central core.
- Adjustable guide and depth gauge.
- Tough Construction for intensive use.
- Swiss manufactured to the highest quality.
- Extraction units with 2 speeds and 4 levels of filtration (including 2 filters to HEPA standards) to 0.3 micron.

**Power:** 100 Watts brushless motor – 24VDC  
**Length:** 200mm  
**Diameter:** 46mm  
**Weight:** 760g  
**Blade Oscillation:** 14,500  
**Noise Level:** 62 dBA  
**Standard:** EC conformity EN 60601  
**Other:** Instant stop button