



Seating

R82 



User height
800–1800 mm



Max. user mass
80 kg



Environment
Indoor & outdoor

Crocodile[®]


The Crocodile is a walking aid for children and young adults to help them into a more upright position. It is placed behind the child and its lightweight design makes it easy to move. The Crocodile is available in 3 sizes.

Highlights about Crocodile



Height adjustable top frame

Easy to adjust the height without using any tools.



Multi-adjustable handles

Positioning the grips close to the body gives better support and a better walking position.



On/off anti-reverse stops

The anti-reverse wheel stops are used for children with less balance and control.



Foldable frame

Makes storage and transporting the Crocodile easy.



Standard

Width and height adjustable handles, lightweight durable aluminium frame - sizes 1 and 2, stable steel and aluminium frame - size 3, upholstery and cover for top frame, solid front and rear wheels, swivel front wheels with directional stabilizer, adjustable drag wheels (rear wheels), on/off anti-reverse wheel stop.

Options

Sizes

- Size 1
- Size 2
- Size 3

Frame colour



Yellow,
Size 1



Orange,
Size 2



Metallic black,
Size 2



Anthracite grey,
Size 3



Crocodile size 1



Crocodile size 3

Measurements

| | | 1 | 2 | 3 |
|--|------------------------------|------------|------------|-------------|
|  | Height from floor to handles | 420-715 mm | 550-840 mm | 750-1045 mm |
|  | Max. user height | 1150 mm | 1500 mm | 1800 mm |
|  | Max. user mass | 30 kg | 45 kg | 80 kg |
|  | Length | 720 mm | 840 mm | 870 mm |
|  | Overall width | 620 mm | 690 mm | 690 mm |
|  | Folded length | 720 mm | 860 mm | 1000 mm |
|  | Folded width | 620 mm | 690 mm | 690 mm |
|  | Folded height | 310 mm | 270 mm | 290 mm |
|  | Minimum turning radius | 1000 mm | 1140 mm | 1140 mm |
|  | Mass device | 5,5 kg | 6,5 kg | 12 kg |

Further measurements at our website

Nice to know

The Crocodile is our hand-held walking aid, allowing children with moderate impairments in standing and walking the possibility to be mobile.

The accessories can be used to provide the exact support that is necessary for the user to feel safe and independent.



Accessories



Hand brakes
86834/86873



Forearm support
86858/59/60



Forearm positioner
869370



Cross bar
86880



Back and side support
8680x



Hip supports
86819




Strap for side supports
Sn4016



Back support
86816



Back support, adjustable
86886



Flip down seat
8682x, 86863



Sling seat
8682x



Handle knobs
86861

For more accessories - visit our website



Bag
86839




Bag
Size 3, 86843



Anti-tips, swing-away
86815 (1-2), 86815-3S (3)



Attendants handle
86866



Watch River walk in the Crocodile size 1. It is also available in size 2 and 3.