

H/S Catheter Set

For enhanced visualization
of uterine pathology

Features & Benefits

- Available with integrated stylet for easy insertion
- Unique placement sheath eliminates the need for a tenaculum
- Comes in both 5 and 7 French
- Ideal for office use
- Latex-free construction
- Made in the USA



Sonohysterography and Hysterosalpingography have been proven to significantly improve the depiction of uterine pathology and determine tubal patency



H/S Procedure Tray

Each Procedure Tray contains:

- Fenestrated drape
- Disposable speculum
- Disposable cervical dilator
- Three swab sticks
- One wide swab
- Sterile wrap
- Two 3" x 3" gauze pads
- Inner tray for Povidone-Iodine solution
- Povidone-Iodine ointment
- Two packs of Povidone-Iodine solution
- 36" Extension tube
- 20 cc syringe with 18-gauge needle
- Lubricating jelly packet
- Sanitary napkin

Complete H/S Procedure Trays


- H/S Procedure Tray - Round 5 French
Product number: 61-5205
- H/S Procedure Tray - Round 7 French
Product number: 61-5207
- H/S Procedure Tray - Round 5 French with Stylet
Product number: 61-3205
- H/S Procedure Tray - Round 7 French with Stylet
Product number: 61-3207

H/S Catheter


The Uterine Infusion Catheter family comes in a standard configuration or with an integrated stylet. Both versions are minimally invasive, easy-to-use devices that ensure a diagnostically enhanced*, yet gentle, patient exam. Designed for use in saline-infusion sonohysterography and hysterosalpingography.

SIS is useful for finding focal lesions in AUB

Pre-Saline Infusion



Post-Saline Infusion



*

Diagnostic Procedure	Sensitivity	Specificity
Transvaginal Ultrasound	75%	76%
Saline Infused Sonohysterography	93%	94%

Kamel, HS, Darwish, AM, Mohamed, SA. Comparison of transvaginal ultrasonography and vaginal sonohysterography in the detection of endometrial polyps. Acta Obstet Gynecol Scand 2000; 79:60

Ordering information

H/S Catheter with Integrated Stylet

Based on the original Ackrad™ round, latex-free, balloon catheter design, this is used when additional catheter rigidity is needed.

- #61-3005 H/S Catheter with integrated Stylet - 5 French, 10 per box
- #61-3007 H/S Catheter with integrated Stylet - 7 French, 10 per box

Standard H/S Catheters

- #61-5005 H/S Catheter, Standard - 5 French, 10 per box
- #61-5007 H/S Catheter, Standard - 7 French, 10 per box

H/Stylet™

- #61-0020 H/Stylet, 10 per box



CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products which aid clinicians in the management and treatment of commonly seen conditions.

203.601.5200 | 800.243.2974 | coopersurgical.com | 95 Corporate Dr. | Trumbull, CT 06611
Also on YouTube at www.YouTube.com/CooperSurgical

All marks bearing the symbol "TM" are trademarks of CooperSurgical, Inc.

© 2011 CooperSurgical, Inc.

Form #80972 • 11/11