

JOBST® CONFIDENCE

Ordering Information

Compression Class / Option

Code

Pieces

JOBST® Confidence

Knee High



Basic style		
CCL 1 (18-21mmHg)	L1-16-04	1
CCL 2 (23-32mmHg)	L2-17-04	1
CCL 3 (34-46mmHg)	L3-18-04	1
Options		
Closed toe	L-A001C	1
SoftFit	L-A016C	1
Slipform	L-A011C	1
Ankle functional zone	L-A010C	1

Options

						
Ankle Functional Zone	Knee Functional Zone	Slipform	Deco Line	Hallux Vallgus (bunion relief zone)*	SoftFit Top Band	Initials
Circumferential functional zone allows greater movement and flexibility at the ankle	Circumferential functional zone allows greater movement and flexibility at the knee	Slipform provides a higher shape to the outside of the garment, at the top, conforming to the limb shape	A decorative seam line can be added to the front of the garment on request	Creates a looser knitting pattern around the ball of the foot to provide relief to painful bunions	Soft and secure top band that has been uniquely knitted into the garment top band for a more comfortable secure fit	Patient initials can be knitted into the garment below the JOBST logo for a personal touch

Compression Class / Option

Code

Pieces

JOBST® Confidence

Thigh High



Basic style		
CCL 1 (18-21mmHg)	L1-16-06	1
CCL 2 (23-32mmHg)	L2-17-06	1
CCL 3 (34-46mmHg)	L3-18-06	1
Options		
Closed toe	L-A001C	1
5cm silicone band	L-A004C	1
Slipform	L-A011C	1
Ankle functional zone	L-A010C	1
Knee functional zone	L-A015C	1

