

JOBST® FarrowWrap®
A-D Below Knee / A-G Thigh High
Ready-to-Wear and Custom Fit Order Form.

TO ORDER PLEASE EMAIL: info@hc21.ie



Please measure all points for a JOBST® FarrowWrap® custom-fit product

E-G

Left:

Right:

D-E

Left:

Right:

A-D

Left:

Right:

Note:

For A-D and D-E, measure posteriorly and follow leg contour.

Note:

For E-G measure medially and follow contours.

G

Groin

(2cm below gluteal fold)

Left:

Right:

F

Upper thigh

Left:

Right:

E

5cm above Patella

Left:

Right:

D

Two finger widths below crease

Left:

Right:

C

Widest calf

Left:

Right:

B²

Midpoint B¹ & C

Left:

Right:

B¹

Base of calf

Left:

Right:

B

Narrowest ankle

Left:

Right:

A¹

Midfoot

Left:

Right:

X

Straight distance

Left:

Right:

A = Floor

Heel to base of little toe

☐ Custom-fit

Brand Compression

Quantity:

LEFT

RIGHT

Lite

Thigh

Classic

Left

Strong

Foot

Remarks

Date:

Purchase Order No.:

Patient Name:

Measured By:

Tel:

Delivery Address:

Email:

Invoice Address:

Please measure all points below for a JOBST® FarrowWrap® ready-to-wear product

	X Small	Small	Medium	Large	X Large	XX Large
Thighpiece (including kneepiece)						
G Groin	60-70 cm	65-75 cm	70-80 cm	75-90 cm	85-100 cm	-
E Above knee	47-52 cm	53-59 cm	60-66 cm	67-73 cm	74-80 cm	-
E-G Short	16-20 cm	16-20 cm	16-20 cm	16-20 cm	16-20 cm	-
E-G Regular	21-26 cm	21-26 cm	21-26 cm	21-26 cm	21-26 cm	-
E-G Tall	27-31 cm	27-31 cm	27-31 cm	27-31 cm	27-31 cm	-
NOTE: for E-G measure medially and follow contours						
Kneepiece (seperate)						
E Above knee	47-52 cm	53-59 cm	60-66 cm	67-73 cm	74-80 cm	81-87 cm
C Mid-Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm	63-73 cm
Legpiece						
C Mid-Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm	-
B Ankle	21-25 cm	25-30 cm	30-36 cm	36-42 cm	42-50 cm	-
A-D Regular	35-39 cm	37-41 cm	39-43 cm	41-45 cm	42-46 cm	-
A-D Tall	40-44 cm	42-46 cm	44-48 cm	46-50 cm	47-51 cm	-
NOTE: for A-D measure posteriorly and follow contours						
Footpiece						
A Mid-foot	22-24 cm	25-27 cm	28-30 cm	31-34 cm	35-40 cm	-
X Regular	15-17 cm	16-19 cm	18-20 cm	19-22 cm	21-24 cm	-
X Tall	18-20 cm	20-22 cm	21-23 cm	23-25 cm	25-27 cm	-
NOTE: for X measure from the heel to the base of the little toe						

Note: 3cm allowance on either side of measurement.

Please use ready-to-wear product code required found on the back of this form

☐ Ready-to-wear

Quantity Left

Product Code

Quantity Right

Product Code

Thigh Piece

Leg Piece

Foot Piece

Remarks

JOBST® FarrowWrap® ready-to-wear direct order form

JOBST® FarrowWrap® Lite (20-30mmHg)

Code	Size	Length
Thighpiece (including kneepiece)		
76655-05	Extra Small	Short
76655-00	Extra Small	Regular
76655-10	Extra Small	Tall
76655-06	Small	Short
76655-01	Small	Regular
76655-11	Small	Tall
76655-07	Medium	Short
76655-02	Medium	Regular
76655-12	Medium	Tall
76655-08	Large	Short
76655-03	Large	Regular
76655-13	Large	Tall
76655-09	Extra Large	Short
76655-04	Extra Large	Regular
76655-14	Extra Large	Tall
Kneepiece (seperate)*		
73367-00	Extra Small	-
73367-01	Small	-
73367-02	Medium	-
73367-03	Large	-
73367-04	Extra Large	-
Legpiece		
76654-00	Extra Small	Regular
76654-05	Extra Small	Tall
76654-01	Small	Regular
76654-06	Small	Tall
76654-02	Medium	Regular
76654-07	Medium	Tall
76654-03	Large	Regular
76654-08	Large	Tall
76654-04	Extra Large	Regular
76654-09	Extra Large	Tall
Footpiece		
76652-05	Extra Small	Regular
76652-00	Extra Small	Long
76652-06	Small	Regular
76652-01	Small	Long
76652-07	Medium	Regular
76652-02	Medium	Long
76652-08	Large	Regular
76652-03	Large	Long
76652-09	Extra Large	Regular
76652-04	Extra Large	Long



JOBST® FarrowWrap® Strong (30-40mmHg)

Code	Size	Length
Thighpiece (including kneepiece)		
76649-05	Extra Small	Short
76649-00	Extra Small	Regular
76649-10	Extra Small	Tall
76649-06	Small	Short
76649-01	Small	Regular
76649-11	Small	Tall
76649-07	Medium	Short
76649-02	Medium	Regular
76649-12	Medium	Tall
76649-08	Large	Short
76649-03	Large	Regular
76649-13	Large	Tall
76649-09	Extra Large	Short
76649-04	Extra Large	Regular
76649-14	Extra Large	Tall
Kneepiece (seperate)*		
73368-00	Extra Small	-
73368-01	Small	-
73368-02	Medium	-
73368-03	Large	-
73368-04	Extra Large	-
73368-05	Extra Extra Large	-
Legpiece		
76648-00	Extra Small	Regular
76648-05	Extra Small	Tall
76648-01	Small	Regular
76648-06	Small	Tall
76648-02	Medium	Regular
76648-07	Medium	Tall
76648-03	Large	Regular
76648-08	Large	Tall
76648-04	Extra Large	Regular
76648-09	Extra Large	Tall
Footpiece		
76646-04	Extra Small	Regular
76646-09	Extra Small	Long
76646-05	Small	Regular
76646-00	Small	Long
76646-06	Medium	Regular
76646-01	Medium	Long
76646-07	Large	Regular
76646-02	Large	Long
76646-08	Extra Large	Regular
76646-03	Extra Large	Long



JOBST® FarrowWrap® Classic (30-40mmHg)

Code	Size	Length
Thighpiece (including kneepiece)		
76658-06	Extra Small	Short
76658-00	Extra Small	Regular
76658-12	Extra Small	Tall
76658-07	Small	Short
76658-01	Small	Regular
76658-13	Small	Tall
76658-08	Medium	Short
76658-02	Medium	Regular
76658-14	Medium	Tall
76658-09	Large	Short
76658-03	Large	Regular
76658-15	Large	Tall
76658-10	Extra Large	Short
76658-04	Extra Large	Regular
76658-16	Extra Large	Tall
Legpiece		
76657-00	Extra Small	Regular
76657-05	Extra Small	Tall
76657-01	Small	Regular
76657-06	Small	Tall
76657-02	Medium	Regular
76657-07	Medium	Tall
76657-03	Large	Regular
76657-08	Large	Tall
76657-04	Extra Large	Regular
76657-09	Extra Large	Tall
Footpiece		
76656-05	Extra Small	Regular
76656-00	Extra Small	Long
76656-06	Small	Regular
76656-01	Small	Long
76656-07	Medium	Regular
76656-02	Medium	Long
76656-08	Large	Regular
76656-03	Large	Long
76656-09	Extra Large	Regular
76656-04	Extra Large	Long

Garment Grip Roll*

Code	Size
Garment Grip Roll	
76672-00	14cm x 80cm
76672-01	14cm x 120cm
76672-02	14cm x 550cm
76672-03	20cm x 120cm
76672-04	30cm x 150cm
76672-05	10cm x 60cm

JOBST® FarrowHybrid®
- supplying compression to the foot and ankle

Code	Size	Colour
ADI (20-30mmHg)		
76664-00	Small (shoe size 2-5)	Beige
76664-01	Medium (shoe size 6-8)	Beige
76664-02	Large (shoe size 9-11)	Beige
76664-03	Extra large (shoe size 12-13)	Beige
76664-12	Small - wide (shoe size 2-5)	Beige
76664-13	Medium - wide (shoe size 6-8)	Beige
76664-14	Large - wide (shoe size 9-11)	Beige
76664-15	Extra large - wide (shoe size 12-13)	Beige
76664-16	Small - wide (shoe size 2-5)	Black
76664-17	Medium - wide (shoe size 6-8)	Black
76664-18	Large - wide (shoe size 9-11)	Black
76664-19	Extra large - wide (shoe size 12-13)	Black
ADII (30-40mmHg)		
76664-04	Small (shoe size 2-5)	Beige
76664-05	Medium (shoe size 6-8)	Beige
76664-06	Large (shoe size 9-11)	Beige
76664-07	Extra large (shoe size 12-13)	Beige
76664-08	Small - wide (shoe size 2-5)	Beige
76664-09	Medium - wide (shoe size 6-8)	Beige
76664-10	Large - wide (shoe size 9-11)	Beige
76664-11	Extra large - wide (shoe size 12-13)	Beige

*Do not apply JOBST® FarrowHybrid if calf circumference exceeds 60cm.
Consider a JOBST® FarrowWrap® footpiece in this circumstance.



A-D Below Knee Liners*

Code	Size
Silver	
76668-00	Small (15-58cm)
76668-02	Large (52-70cm)
TG Soft	
76669-00	Small (<40cm)
76669-02	Medium (40-70cm)
76669-04	Large (70-120cm)

A-G Thigh High Liners*

Code	Size
TG Soft	
76670-02	Medium (40-70cm)
76670-04	Large (70-120cm)

Foam Squares Pack*

Code	Size
Foam Squares Pack	
76675-01	Small (8x8cm)
76675-02	Large (10x10cm)

NB The thighpiece comes with kneepiece included, however, the kneepiece can be ordered separately.

Date: _____ Purchase Order No.: _____ Patient Name: _____

Measured By: _____ Tel: _____ Email: _____

Delivery Address: _____ Invoice Address: _____

