


JOBST® FarrowWrap® Lite Armiece and Gauntlet

Ready-to-wear order form



TO ORDER PLEASE EMAIL: info@hc21.ie



Axilla Left: Right:

Circumference and length measurements in centimetres (cms)

Elbow Left: Right:

Left: Right:

Wrist Left: Right:

Hand A: B:

A - Palm at base of fingers
B - Palm at fold of the thumb

NOTE: Measure length from the outside of the wrist to the axilla (G) with a slight bend in arm

Date: _____ Purchase Order No.: _____ Patient Name: _____

Measured By: _____ Tel: _____ Email: _____

Delivery Address: _____ Invoice Address: _____

Code	Option	Size	Length	Colour
Armiece				
76650-11	Left	Small	Short	Beige
76650-05	Left	Small	Regular	Beige
76650-35	Left	Small	Long	Beige
76650-13	Left	Medium	Short	Beige
76650-07	Left	Medium	Regular	Beige
76650-01	Left	Medium	Long	Beige
76650-15	Left	Large	Short	Beige
76650-09	Left	Large	Regular	Beige
76650-03	Left	Large	Long	Beige
76650-29	Left	Small	Short	Black
76650-23	Left	Small	Regular	Black
76650-13	Left	Small	Long	Black
76650-31	Left	Medium	Short	Black
76650-25	Left	Medium	Regular	Black
76650-19	Left	Medium	Long	Black
76650-33	Left	Large	Short	Black
76650-27	Left	Large	Regular	Black
76650-21	Left	Large	Long	Black
76650-12	Right	Small	Short	Beige
76650-06	Right	Small	Regular	Beige
76650-00	Right	Small	Long	Beige
76650-14	Right	Medium	Short	Beige
76650-08	Right	Medium	Regular	Beige
76650-02	Right	Medium	Long	Beige
76650-16	Right	Large	Short	Beige
76650-10	Right	Large	Regular	Beige
76650-04	Right	Large	Long	Beige
76650-30	Right	Small	Short	Black
76650-24	Right	Small	Regular	Black
76650-18	Right	Small	Long	Black
76650-32	Right	Medium	Short	Black
76650-26	Right	Medium	Regular	Black
76650-20	Right	Medium	Long	Black
76650-34	Right	Large	Short	Black
76650-28	Right	Large	Regular	Black
76650-22	Right	Large	Long	Black
Gauntlet				
76653-00	Ambidextrous	Extra small	-	Beige
76653-01	Ambidextrous	Small	-	Beige
76653-02	Ambidextrous	Medium	-	Beige
76653-03	Ambidextrous	Large	-	Beige

Armiece

Size: Small Medium Large **Side:** Left Right

Length: Short Regular Long **Colour:** Tan Black

Compression: Lite (20-30 mmHg)

Code: -

Quantity required:

Repeat prescription required every months

Gauntlet

Size: X Small Small Medium Large **Colour:** Tan

Compression: Lite (20-30 mmHg)

Code: -

Quantity required:

Repeat prescription required every months



	Wrist	Elbow	Axilla		Length
Small	14cm - 18cm	20cm - 27cm	22cm - 31cm	Short	40cm - 43cm
Medium	16cm - 21cm	25cm - 34cm	29cm - 39cm	Regular	44cm - 47cm
Large	19cm - 25cm	30cm - 40cm	32cm - 45cm	Long	48cm - 51cm

Gauntlet

	X Small	Small	Medium	Large
A - Palm at base of fingers	< 18cm	18-20cm	20-22cm	> 22cm
B - Palm at fold of thumb	< 19cm	19-21cm	21-23cm	> 24cm

