


JOBST® FarrowWrap®
A-G Thigh High
Ready-to-wear order form

TO ORDER PLEASE EMAIL: info@hc21.ie



A pair of appropriately sized liners will be provided free of charge with each JOBST® FarrowWrap® lower limb piece.
Below knee comes with silver liners, thigh high comes with TG Soft liners.

Customer Service: **1890 777 444** Compression Therapy Helpline: **0800 389 8424**
Email: info@hc21.ie Website: www.hc21.eu

Date: _____ Purchase Order No.: _____

Measured By: _____

Delivery Address: _____

Patient Name: _____

Tel: _____ Email: _____

Invoice Address: _____

	X Small	Small	Medium	Large	X Large
Thighpiece (including kneepiece)					
G Groin	60-70 cm	65-75 cm	70-80 cm	75-90 cm	85-100 cm
E Above knee	47-52 cm	53-59 cm	60-66 cm	67-73 cm	74-80 cm
E-G Short	16-20 cm	16-20 cm	16-20 cm	16-20 cm	16-20 cm
E-G Regular	21-26 cm	21-26 cm	21-26 cm	21-26 cm	21-26 cm
E-G Tall	27-31 cm	27-31 cm	27-31 cm	27-31 cm	27-31 cm
NOTE: for E-G measure medially and follow contours					
Legpiece					
C Mid-Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm
B Ankle	21-25 cm	25-30 cm	30-36 cm	36-42 cm	42-50 cm
A-D Regular	35-39 cm	37-41 cm	39-43 cm	41-45 cm	42-46 cm
A-D Tall	40-44 cm	42-46 cm	44-48 cm	46-50 cm	47-51 cm
NOTE: for A-D measure posteriorly and follow contours					
Footpiece					
A¹ Mid-foot	22-24 cm	25-27 cm	28-30 cm	31-34 cm	35-40 cm
X Regular	15-17 cm	16-19 cm	18-20 cm	19-22 cm	21-24 cm
X Tall	18-20 cm	20-22 cm	21-23 cm	23-25 cm	25-27 cm
NOTE: for X measure from the heel to the base of the little toe					

Note: 3cm allowance on either side of measurement.

☐ **Thighpiece**

Size: ☐ X Small ☐ Small ☐ Medium ☐ Large ☐ X Large

Length E-G: ☐ Short ☐ Regular ☐ Tall

Fabric: ☐ Lite ☐ Classic ☐ Strong

Code: -

Quantity Required:

Repeat prescription required every months

☐ **Legpiece**

Size: ☐ X Small ☐ Small ☐ Medium ☐ Large ☐ X Large

Length A-D: ☐ Regular ☐ Tall

Fabric: ☐ Lite ☐ Classic ☐ Strong

Code: -

Quantity Required:

Repeat prescription required every months

☐ **Footpiece**

Size: ☐ X Small ☐ Small ☐ Medium ☐ Large ☐ X Large

Length X: ☐ Regular ☐ Long

Fabric: ☐ Lite ☐ Classic ☐ Strong

Code: -

Quantity Required:

Repeat prescription required every months

☐ **JOBST® FarrowHybrid®**

Size: ☐ Small ☐ Medium ☐ Large ☐ X Large

Width: ☐ Regular* ☐ Wide**

Compression: 20-30mmHg - ☐ Taupe ☐ Black*** 30-40mmHg - ☐ Sand

Code: -

Quantity Required:

Repeat prescription required every months

*Regular: 20-27 cm **A¹** **Wide: 28-35 cm **A¹** ***Available in wide variant only

JOBST® FarrowWrap® Order Codes

JOBST® FarrowWrap® Lite (20-30mmHg)

Code	Size	Length	Colour	PIP Code
Thighpiece (including kneepiece)				
76655-05	Extra Small	Short	Beige	
76655-00	Extra Small	Regular	Beige	
76655-10	Extra Small	Tall	Beige	
76655-06	Small	Short	Beige	
76655-01	Small	Regular	Beige	
76655-11	Small	Tall	Beige	
76655-07	Medium	Short	Beige	
76655-02	Medium	Regular	Beige	
76655-12	Medium	Tall	Beige	
76655-08	Large	Short	Beige	
76655-03	Large	Regular	Beige	
76655-13	Large	Tall	Beige	
76655-09	Extra Large	Short	Beige	
76655-04	Extra Large	Regular	Beige	
76655-14	Extra Large	Tall	Beige	
Legpiece				
76654-00	Extra Small	Regular	Beige	
76654-05	Extra Small	Tall	Beige	
76654-01	Small	Regular	Beige	
76654-06	Small	Tall	Beige	
76654-02	Medium	Regular	Beige	
76654-07	Medium	Tall	Beige	
76654-03	Large	Regular	Beige	
76654-08	Large	Tall	Beige	
76654-04	Extra Large	Regular	Beige	
76654-09	Extra Large	Tall	Beige	
Footpiece				
76652-05	Extra Small	Regular	Beige	
76652-00	Extra Small	Long	Beige	
76652-06	Small	Regular	Beige	
76652-01	Small	Long	Beige	
76652-07	Medium	Regular	Beige	
76652-02	Medium	Long	Beige	
76652-08	Large	Regular	Beige	
76652-03	Large	Long	Beige	
76652-09	Extra Large	Regular	Beige	
76652-04	Extra Large	Long	Beige	



JOBST® FarrowWrap® Classic (30-40mmHg)

Code	Size	Length	Colour	PIP Code
Thighpiece (including kneepiece)				
76658-06	Extra Small	Short	Beige	
76658-00	Extra Small	Regular	Beige	
76658-12	Extra Small	Tall	Beige	
76658-07	Small	Short	Beige	
76658-01	Small	Regular	Beige	
76658-13	Small	Tall	Beige	
76658-08	Medium	Short	Beige	
76658-02	Medium	Regular	Beige	
76658-14	Medium	Tall	Beige	
76658-09	Large	Short	Beige	
76658-03	Large	Regular	Beige	
76658-15	Large	Tall	Beige	
76658-10	Extra Large	Short	Beige	
76658-04	Extra Large	Regular	Beige	
76658-16	Extra Large	Tall	Beige	
Legpiece				
76657-00	Extra Small	Regular	Beige	
76657-05	Extra Small	Tall	Beige	
76657-01	Small	Regular	Beige	
76657-06	Small	Tall	Beige	
76657-02	Medium	Regular	Beige	
76657-07	Medium	Tall	Beige	
76657-03	Large	Regular	Beige	
76657-08	Large	Tall	Beige	
76657-04	Extra Large	Regular	Beige	
76657-09	Extra Large	Tall	Beige	
Footpiece				
76656-05	Extra Small	Regular	Beige	
76656-00	Extra Small	Long	Beige	
76656-06	Small	Regular	Beige	
76656-01	Small	Long	Beige	
76656-07	Medium	Regular	Beige	
76656-02	Medium	Long	Beige	
76656-08	Large	Regular	Beige	
76656-03	Large	Long	Beige	
76656-09	Extra Large	Regular	Beige	
76656-04	Extra Large	Long	Beige	




JOBST® FarrowWrap® Strong (30-40mmHg)

Code	Size	Length	Colour	PIP Code
Thighpiece (including kneepiece)				
76649-05	Extra Small	Short	Beige	
76649-00	Extra Small	Regular	Beige	
76649-10	Extra Small	Tall	Beige	
76649-06	Small	Short	Beige	
76649-01	Small	Regular	Beige	
76649-11	Small	Tall	Beige	
76649-07	Medium	Short	Beige	
76649-02	Medium	Regular	Beige	
76649-12	Medium	Tall	Beige	
76649-08	Large	Short	Beige	
76649-03	Large	Regular	Beige	
76649-13	Large	Tall	Beige	
76649-09	Extra Large	Short	Beige	
76649-04	Extra Large	Regular	Beige	
76649-14	Extra Large	Tall	Beige	
Legpiece				
76648-00	Extra Small	Regular	Beige	
76648-05	Extra Small	Tall	Beige	
76648-01	Small	Regular	Beige	
76648-06	Small	Tall	Beige	
76648-02	Medium	Regular	Beige	
76648-07	Medium	Tall	Beige	
76648-03	Large	Regular	Beige	
76648-08	Large	Tall	Beige	
76648-04	Extra Large	Regular	Beige	
76648-09	Extra Large	Tall	Beige	
Footpiece				
76646-04	Extra Small	Regular	Beige	
76646-09	Extra Small	Long	Beige	
76646-05	Small	Regular	Beige	
76646-00	Small	Long	Beige	
76646-06	Medium	Regular	Beige	
76646-01	Medium	Long	Beige	
76646-07	Large	Regular	Beige	
76646-02	Large	Long	Beige	
76646-08	Extra Large	Regular	Beige	
76646-03	Extra Large	Long	Beige	

JOBST® FarrowHybrid®* - supplying compression to the foot and ankle

Code	Size	Colour	PIP Code	Code	Size	Colour	PIP Code
ADI (20-30mmHg)				ADII (30-40mmHg)			
76664-00	Small (shoe size 2-5)	Beige		76664-04	Small (shoe size 2-5)	Beige	
76664-01	Medium (shoe size 6-8)	Beige		76664-05	Medium (shoe size 6-8)	Beige	
76664-02	Large (shoe size 9-11)	Beige		76664-06	Large (shoe size 9-11)	Beige	
76664-03	Extra large (shoe size 12-13)	Beige		76664-07	Extra large (shoe size 12-13)	Beige	
76664-12	Small - wide (shoe size 2-5)	Beige		76664-08	Small - wide (shoe size 2-5)	Beige	
76664-13	Medium - wide (shoe size 6-8)	Beige		76664-09	Medium - wide (shoe size 6-8)	Beige	
76664-14	Large - wide (shoe size 9-11)	Beige		76664-10	Large - wide (shoe size 9-11)	Beige	
76664-15	Extra large - wide (shoe size 12-13)	Beige		76664-11	Extra large - wide (shoe size 12-13)	Beige	
76664-16	Small - wide (shoe size 2-5)	Black		*Do not apply JOBST® FarrowHybrid® if calf circumference exceeds 60cm. Consider a JOBST® FarrowWrap® footpiece in this circumstance.			
76664-17	Medium - wide (shoe size 6-8)	Black					
76664-18	Large - wide (shoe size 9-11)	Black					
76664-19	Extra large - wide (shoe size 12-13)	Black					



Measuring instructions

Thighpiece - 3 simple measuring points

1. Measure 5cm above the patella **E** and the groin **G** circumference (as shown on diagram)
2. Use chart below to determine size at the top of the column (if 2 measurements do not fit in same column the smaller size should accommodate reducing oedema)
3. Measure the length of the thigh medially between **E** and **G** following the contour of the leg. Use the sizing chart to select short, regular or tall to determine size

Legpiece - 3 simple measuring points

1. Measure the widest calf **C** and narrowest ankle **B** circumference (as shown on diagram)
2. Use chart below to determine size at the top of the column (if 2 measurements do not fit in same column the smaller size should accommodate reducing the oedema)
3. Measure the length of the leg posteriorly from the floor to **D** measurement (2cm below the crease of the knee). Follow the contour of the leg. Use the sizing chart to select regular or tall to determine size

Footpiece - 2 simple measuring points

1. Measure **A** circumference (as shown on diagram)
2. Use chart below to determine size at the top of the column
3. Measure the **X** length of the foot. Use the sizing chart to select regular or long to determine size