



Mityvac® and Mystic® II

Softer, safer vacuum assisted delivery systems

Enhancing a winner —
the Mityvac line of vacuum
assisted delivery systems

Features

■ Gentle

MitySoft Bell cup protects vaginal wall and cradles baby's head

■ Low Profile

Flexible Mystic M-Style Mushroom® cup flexes up to 90° in any direction to facilitate use in non occiput anterior positions

■ Ease of Use

Prominent trigger release disengages vacuum

■ Safety

Handle rotates without transferring to the cup

■ Self Limiting

Mystic II's is specially designed vacuum will not exceed 5 cmHg (580 mm Hg)

Advanced Mystic II pump and cups


M-Style® low-profile
Mushroom® cup

MitySoft® Bell cup

Extended shaft

Vacuum release

Auxiliary release/
connection valve



CooperSurgical

*Keeping You at the Forefront
of Women's Health Care™*

ORDERING INFORMATION

Product #	Description of Cups and Accessories	Diameter	Order Quantity
10004	Pearl Edge® Bell Cup with 4" tubing and filter	60 mm	12 per box
10007LP	M-Style® Mushroom® Cup with 4" tubing and filter	50 mm	12 per box
10008	Super M-Style® Mushroom® Cup with 4" tubing and filter	56 mm	12 per box
10020	MitySoft® Bell Cup with 4" tubing and filter	64 mm	12 per box
10137	M-Select® Mushroom® Cup with 4" tubing and filter	50 mm	12 per box
10500	Reusable Vacuum Delivery Cup	70 mm	1 per box
Pumps			
10022	Mityvac® Reusable Pump	N/A	1 per box
10057	Mystic® II Pump with M-Style® Mushroom® Cup	50 mm	12 per box
10058	Mystic® II Pump with MitySoft® Bell Cup	64 mm	12 per box
Accessories			
16177	Mityvac® Reusable Pump Replacement Gauge	N/A	1 per box
16710	Mityvac® Reusable Pump Maintenance Kit	N/A	1 per box



www.coopersurgical.com

To find out more about the Mityvac® and Mystic® II, or to place an order call 800.243.2974.

The marks bearing the symbol "®" are registered trademarks of CooperSurgical.

Mityvac® and Mystic® II

95 Corporate Drive, Trumbull, CT 06611
800.243.2974 • 203.601.5200