



AirLife™
Noninvasive ventilation masks

AirLife™ noninvasive ventilation masks


Our noninvasive ventilation (NIV) full face and nasal masks offer clinicians a reliable option while providing patients comfort and security.

The adjustable four-point headgear enables ease in fit and adjustment while the flexible swivel elbow reduces torque applied to the seal between the mask and patient, helping to minimize pressure points.

Our masks feature a lightweight silicone design that contours to the patient's face for a comfortable fit, and the quick-release clip eliminates the need to refit after removal.

The odorless, easy-to-clean masks include both vented and non-vented options and are available in three sizes to fit a wide range of patients.

The mask flexes when the circuit is pulled.



Mask when flexed



Mask not flexed



Full face mask

Key features

- Flexible swivel connector enables ease in positioning of the circuit
- Comfortable four-point headgear with crown strap for better stability and sealing
- Adjustable forehead rest for comfort and stability
- Soft seal that conforms to the face, reducing need to overtighten
- Quick-release strap
- Made without DEHP

AirLife NIV masks (10/cs.)

Cat no.	Product description
NIV040S	Full face mask, vented with anti-asphyxia valve (AAV), small
NIV040M	Full face mask, vented with AAV, medium
NIV040L	Full face mask, vented with AAV, large
NIV043S	Full face mask, non-vented with AAV, small
NIV043M	Full face mask, non-vented with AAV, medium
NIV043L	Full face mask, non-vented with AAV, large
NIV041S	Full face mask, non-vented, small
NIV041M	Full face mask, non-vented, medium
NIV041L	Full face mask, non-vented, large
NIV042S	Nasal mask, vented, small
NIV042M	Nasal mask, vented, medium
NIV042L	Nasal mask, vented, large
NIV044S	Nasal mask, non-vented with AAV, small
NIV044M	Nasal mask, non-vented with AAV, medium
NIV044L	Nasal mask, non-vented with AAV, large

For more information, please call Customer Service at **800.323.9088**.



Aquilant

A HC21 Group Company

Unit B1-B2, Bond Close, Basingstoke, RG24 8PZ

E: contactus@hc21.eu

T: +44 (0) 1256 306506

www.healthcare21.eu

GLOBAL HEADQUARTERS

Vyaire Medical, Inc.

26125 North Riverwoods Blvd

Mettawa, IL 60045, USA

vyaire.com

Rx ONLY For U.S. distribution only.

© 2020 Vyaire. Vyaire, the Vyaire logo and AirLife are trademarks or registered trademarks of Vyaire Medical, Inc., or one of its affiliates. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-US-2000139

vyaire[™]
M E D I C A L