




CE
0120


Multi electrode for large & Multiple lesions

Octopus RF ELECTRODE



Octopus RF electrodes are designed to activate three active tips with 200 watts simultaneously for creating a large ablation zone.

For radiofrequency ablation of multiple lesions, Octopus RF electrodes are able to be separated to three parts to create ablation zone respectively.



Octopus RF Electrode SPECIFICATION

(Unit: mm)

Ordering No.	Diameter (Gauge)	Length	Tip exposure length
17-15P20	17G	150	20
17-15P25	17G	150	25
17-15P30	17G	150	30
17-20P20	17G	200	20
17-20P25	17G	200	25
17-20P30	17G	200	30
15-15P25	15G	150	25
15-15P30	15G	150	30
15-20P25	15G	200	25
15-20P30	15G	200	30