


# Single Use Endoscope Cleaning Brush Set

A high quality endoscope cleaning brush set designed to ensure maximum effectiveness and efficiency

Convenient pack containing:

1 x channel cleaning brush

1 x Port & Valve cleaning brush



**Helping you to deliver  
world class endoscopy services**

Designed and  
manufactured in the UK

# Channel Cleaning Brush Set Features:

- Single use
- 3 different sizes available to ensure maximum cleaning effectiveness
- Colour coded to simplify correct brush size selection
- Innovative swage design to reduce the wear and tear on the endoscope channel
- High quality nylon bristle construction, ensures effective debris removal without causing channel damage
- Port & Valve housing stubby brush with comfortable grooved, slip resistant design
- Two brush pack to simplify ordering and stock management
- 2 x traceability stickers to comply with disinfection guidelines
- High quality components and UK based manufacturing reducing time consuming and costly product faults



Description	Product Code	Quantity	Sheath Diameter	Bristle Diameter	Channel Size	Length
Channel Cleaning Brush Set (Red)	AQCCB/2.8-2500	50	2mm	5mm	2.8mm – 4.8mm	2500mm
Channel Cleaning Brush Set (Teal)	AQCCB/1.8-2500	50	2mm	5mm	1.8mm – 2.4mm	2500mm
Channel Cleaning Brush Set (Blue)	AQCCB/1.2-1000	50	1.2mm	5mm	1.2mm – 1.5mm	1000mm
Port & Valve Housing Stubby Brush	AQPVB/2-11	50	N/A	Large 5mm – 11mm Small 2mm – 5mm	N/A	140mm

**Customer Order Line:** T + 44 (0) 1256 365 456    [contactus@aquilantendoscopy.com](mailto:contactus@aquilantendoscopy.com)

## Aquilant Specialist Healthcare Services

Aquilant House, Unit B1-B2, Bond Close, Kingsland Business Park,  
Basingstoke, Hampshire, RG24 8PZ, United Kingdom

T + 44 (0) 1256 365 456    [contactus@aquilantendoscopy.com](mailto:contactus@aquilantendoscopy.com)    [www.aquilantendoscopy.com](http://www.aquilantendoscopy.com)