



# Esophageal Dilatation Systems

# M-Flex® Blue Silicone Bougies



There’s nothing more flexible than an M-Flex.®

Our silicone bougies feature added flexibility that some users prefer. Tungsten-filled and not made with natural rubber latex, there’s nothing quite like the M-Flex.®

## Benefits

- Ultra-flexible tapered tip (Maloney)
- Completely radiopaque from distal tip to proximal end
- American and European depth markings
- Mercury-free and not made with natural rubber latex
- Autoclavable for quick reprocessing

## Features

- 21 sizes: 20 FR to 60 FR (in 2 FR increments)
- Choose a full set, mini set or individual pieces
- Three-year warranty on product and workmanship
- Easy disposal
- No charge trade-in program for proper disposal of expired mercury-filled bougies (with purchase of new product)



Store and transport M-Flex® bougies in our storage case which holds up to two full sets.

## M-Flex® Dilatation Systems

| Set                                 | Part No. | Contents                                                  |
|-------------------------------------|----------|-----------------------------------------------------------|
| Maloney (tapered) 21-piece full set | 1206-00  | 21 bougies, 20–60 FR                                      |
| Maloney (tapered) 10-piece mini set | 1206-10  | 10 bougies, 36–54 FR                                      |
| Hurst (blunt) 21-piece full set     | 1207-00  | 21 bougies, 20–60 FR                                      |
| Hurst (blunt) 10-piece mini set     | 1207-10  | 10 bougies, 36–54 FR                                      |
| Bougie Storage Case                 | 1210-01  | Heavy-duty storage for up to 46 bougies. 32"L x 7"W x 4"H |




Our storage case allows easy identification of bougie size. Removing plexiglass window provides easy access.



End caps display French size.



Extra flexible tip.



Top: Maloney tip. Bottom: Hurst tip.

## M-Flex® Blue Silicone Bougies

| FR (French) Size | (mm) | Maloney (tapered) Tip Part No. | Hurst (blunt) Tip Part No. |
|------------------|------|--------------------------------|----------------------------|
| 20               | 6.6  | 1206-20                        | 1207-20                    |
| 22               | 7.4  | 1206-22                        | 1207-22                    |
| 24               | 8    | 1206-24                        | 1207-24                    |
| 26               | 8.6  | 1206-26                        | 1207-26                    |
| 28               | 9.3  | 1206-28                        | 1207-28                    |
| 30               | 10   | 1206-30                        | 1207-30                    |
| 32               | 10.7 | 1206-32                        | 1207-32                    |
| 34               | 11.2 | 1206-34                        | 1207-34                    |
| 36               | 12   | 1206-36                        | 1207-36                    |
| 38               | 12.5 | 1206-38                        | 1207-38                    |
| 40               | 13   | 1206-40                        | 1207-40                    |
| 42               | 14   | 1206-42                        | 1207-42                    |
| 44               | 14.5 | 1206-44                        | 1207-44                    |
| 46               | 15.2 | 1206-46                        | 1207-46                    |
| 48               | 16   | 1206-48                        | 1207-48                    |
| 50               | 16.5 | 1206-50                        | 1207-50                    |
| 52               | 17.2 | 1206-52                        | 1207-52                    |
| 54               | 18   | 1206-54                        | 1207-54                    |
| 56               | 18.5 | 1206-56                        | 1207-56                    |
| 58               | 19.2 | 1206-58                        | 1207-58                    |
| 60               | 20   | 1206-60                        | 1207-60                    |

Hurst Bougie

Maloney Bougie

American Depth Markings

European Depth Markings

# WeightRight® Mercury-Free Bougies

## Heavy-weight control and convenience combined.

Our tungsten-filled bougie mimics the weight and feel of mercury-filled bougies without the concerns of mercury usage and disposal. It's the weight you expect, and the convenience you demand.

### Benefits

- Completely radiopaque from distal tip to proximal end
- Flexible material does not attract particulate
- American and European depth markings
- Mercury-free and not made with natural rubber latex

### Features

- 23 sizes: 16 FR to 60 FR (in 2 FR increments)
- Choose a full set, mini set or individual pieces
- Three-year warranty on product and workmanship
- Easy disposal
- No charge trade-in program for proper disposal of expired mercury-filled bougies (with purchase of new product)

### WeightRight® Dilatation Systems

| Set                                 | Part No. | Contents                                                  |
|-------------------------------------|----------|-----------------------------------------------------------|
| Maloney (tapered) 23-piece full set | 1208-00  | 23 bougies, 16-60 FR                                      |
| Maloney (tapered) 10-piece mini set | 1208-10  | 10 bougies, 36-54 FR                                      |
| Hurst (blunt) 23-piece full set     | 1209-00  | 23 bougies, 16-60 FR                                      |
| Hurst (blunt) 10-piece mini set     | 1209-10  | 10 bougies, 36-54 FR                                      |
| Bougie Storage Case                 | 1210-01  | Heavy-duty storage for up to 46 bougies. 32"L x 7"W x 4"H |



Store and transport WeightRight® bougies in our storage case.



Our storage case allows easy identification of bougie size. Removing plexiglass window provides easy access.



Hang bougie vertically using the hole in the end cap.



Top: Maloney tip. Bottom: Hurst tip.

## WeightRight® Mercury-free Bougies

| FR (French) Size | (mm) | Maloney (tapered) Tip Part No. | Hurst (blunt) Tip Part No. |
|------------------|------|--------------------------------|----------------------------|
| 16               | 5.3  | 1208-16                        | 1209-16                    |
| 18               | 6    | 1208-18                        | 1209-18                    |
| 20               | 6.6  | 1208-20                        | 1209-20                    |
| 22               | 7.4  | 1208-22                        | 1209-22                    |
| 24               | 8    | 1208-24                        | 1209-24                    |
| 26               | 8.6  | 1208-26                        | 1209-26                    |
| 28               | 9.3  | 1208-28                        | 1209-28                    |
| 30               | 10   | 1208-30                        | 1209-30                    |
| 32               | 10.7 | 1208-32                        | 1209-32                    |
| 34               | 11.2 | 1208-34                        | 1209-34                    |
| 36               | 12   | 1208-36                        | 1209-36                    |
| 38               | 12.5 | 1208-38                        | 1209-38                    |
| 40               | 13   | 1208-40                        | 1209-40                    |
| 42               | 14   | 1208-42                        | 1209-42                    |
| 44               | 14.5 | 1208-44                        | 1209-44                    |
| 46               | 15.2 | 1208-46                        | 1209-46                    |
| 48               | 16   | 1208-48                        | 1209-48                    |
| 50               | 16.5 | 1208-50                        | 1209-50                    |
| 52               | 17.2 | 1208-52                        | 1209-52                    |
| 54               | 18   | 1208-54                        | 1209-54                    |
| 56               | 18.5 | 1208-56                        | 1209-56                    |
| 58               | 19.2 | 1208-58                        | 1209-58                    |
| 60               | 20   | 1208-60                        | 1209-60                    |

Hurst Bougie

Maloney Bougie

American Depth Markings

European Depth Markings

# SafeGuide® Over the Guidewire Esophageal Dilators

## Flexibility and control for exacting procedures.

At times, there is need for greater accuracy. That's why there's SafeGuide®, a remarkably flexible over the guidewire solution that ensures exact placement of the dilator. Position the guidewire with the assistance of an endoscope and direct the SafeGuide® dilator where you need it.

### Benefits

- Completely radiopaque from distal tip to proximal end
- Flexible material does not attract particulate
- American and European depth markings
- Not made with natural rubber latex
- Dilators are reusable

### Features

- 16 sizes: 15 FR to 60 FR (in 3 FR increments)
- Choose a full set, mini set or individual pieces
- Easy to clean with optional brushes and adapters
- With proper care, reusable SafeGuide® dilators provide years of service


### SafeGuide® Dilatation Systems

| Set               | Part No.                          | Contents                                                                                                                                              |
|-------------------|-----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| 16-piece full set | 1214-00                           | 16 dilators, 15-60 FR, 1 storage case, 1 marked spring tip guidewire, 1 reusable cleaning brush, 25 disposable cleaning adapters                      |
| 7-piece mini set  | 1214-10                           | 7 dilators (15, 21, 27, 33, 39, 42, 45 FR), 1 storage case, 1 marked spring tip guidewire, 1 reusable cleaning brush, 25 disposable cleaning adapters |
| Custom mini set   | Contact your sales representative | 7 dilators of your choice, 1 storage case, 1 marked spring tip guidewire, 1 reusable cleaning brush, 25 disposable cleaning adapters                  |

Store and transport SafeGuide® dilators in a convenient storage case. A labeled slot holds each dilator secure. The case holds a full set of 16 dilators.



The single-piece construction of the SafeGuide® Spring Tip Guidewire resists breakage and flexes easily.



### SafeGuide® Over the Guidewire Esophageal Dilators

| FR (French) Size | (mm) | Part No. |
|------------------|------|----------|
| 15               | 5    | 1214-15  |
| 18               | 6    | 1214-18  |
| 21               | 7    | 1214-21  |
| 24               | 8    | 1214-24  |
| 27               | 9    | 1214-27  |
| 30               | 10   | 1214-30  |
| 33               | 11   | 1214-33  |
| 36               | 12   | 1214-36  |
| 39               | 13   | 1214-39  |
| 42               | 14   | 1214-42  |
| 45               | 15   | 1214-45  |
| 48               | 16   | 1214-48  |
| 51               | 17   | 1214-51  |
| 54               | 18   | 1214-54  |
| 57               | 19   | 1214-57  |
| 60               | 20   | 1214-60  |

### Accessories



Storage case indicates size.

### SafeGuide® Accessories

| Set                          | Part No. | Description                                           | Qty/Case |
|------------------------------|----------|-------------------------------------------------------|----------|
| SafeGuide® Storage Case      | 1214-01  | Holds up to 16 SafeGuide® Dilators. 36"L x 10"W x 4"H | 1        |
| Marked Spring Tip Guidewires | 1214-02  | Guides the dilator. 210cm                             | 2        |
| Disposable Cleaning Brushes  | 1214-03D | Cleans the inner lumen of the dilator                 | 3        |
| Reusable Cleaning Brushes    | 1214-03R | Cleans the inner lumen of the dilator                 | 3        |
| Disposable Cleaning Adapters | 1214-04  | Used to flush the inner lumen of the dilator          | 25       |



Clearly marked reusable and single-use cleaning brushes.



Insert cleaning brush into the distal end and push through to the other side for complete inner lumen cleaning.



Disposable cleaning adapter.

## Is it time to replace your esophageal bougies or wire-guided dilators?

### Replace bougies when you notice:

- End cap separating from bougie
- Cuts, cracks or nicks
- Flaking or peeling of the outer coating
- Faded or missing markings
- Expired or missing expiration dates

### Replace wire-guided dilators when you notice:

- Loss of flexibility
- Cuts, cracks or nicks
- Stretching in the inner lumen
- Faded or missing markings
- Curved dilator or curled tapered tip
- Expired



**Esophageal  
Dilators**



**CleanFreak®  
Cleaning  
Products**



**Endoscope  
Care &  
Accessories**

At Medovations, we are committed to providing you with high-quality, durable products for your gastrointestinal procedures. See the complete line of GI solutions at **Medovations.com**

Learn more from your personal product representative at **800-558-6408** or **Medovations.com**



**medovations®**

Medovations, Inc. | 102 East Keefe Avenue | Milwaukee, Wisconsin 53212 USA  
800.558.6408 or 414.265.7620 | Fax 414.265.7628  
Orders: **orders@medovations.com** | General Inquiries: medo@medovations.com